

Endomech and Energy

Product Catalogue

Nordics

Our Credo

We believe our first responsibility is to the patients, doctors and nurses, to mothers and fathers and all others who use our products and services. In meeting their needs everything we do must be of high quality. We must constantly strive to provide value, reduce our costs and maintain reasonable prices. Customers' orders must be serviced promptly and accurately. Our business partners must have an opportunity to make a fair profit.

We are responsible to our employees who work with us throughout the world. We must provide an inclusive work environment where each person must be considered as an individual. We must respect their diversity and dignity and recognize their merit. They must have a sense of security, fulfillment and purpose in their jobs. Compensation must be fair and adequate, and working conditions clean, orderly and safe. We must support the health and well-being of our employees and help them fulfill their family and other personal responsibilities. Employees must feel free to make suggestions and complaints. There must be equal opportunity for employment, development and advancement for those qualified. We must provide highly capable leaders and their actions must be just and ethical.

We are responsible to the communities in which we live and work and to the world community as well. We must help people be healthier by supporting better access and care in more places around the world. We must be good citizens – support good works and charities, better health and education, and bear our fair share of taxes. We must maintain in good order the property we are privileged to use, protecting the environment and natural resources.

Our final responsibility is to our stockholders. Business must make a sound profit. We must experiment with new ideas. Research must be carried on, innovative programs developed, investments made for the future and mistakes paid for. New equipment must be purchased, new facilities provided and new products launched. Reserves must be created to provide for adverse times. When we operate according to these principles, the stockholders should realize a fair return.

Johnson & Johnson



Contact us

Sweden:

Tel: 08-626 22 03
E-post: cservse@its.jnj.com

Denmark:

Tel: 45 94 82 04
E-post: cservdk@its.jnj.com

Finland:

Puh: 0207-531 502
s-posti: cservfi@its.jnj.com

Norway:

Tel: 67 51 86 00
E-post: ordre@ortomedic.no

Customer Connect

J&J Customer Connect is a global portal for the Johnson & Johnson Family of Companies that lets you place orders, obtain order status, access our product catalog and more.


Our enhanced site offers:

Easier navigation

Easier-to-find order details


Mobile device optimization

J&J Customer Connect portal offers a 24/7 direct information access of daily inquiries and generates time-saving in the daily follow-up.



[HTTPS://JJCUSTOMERCONNECT.COM](https://jjcustomerconnect.com)

Our goal is to provide you with an easy-to-use portal to place orders and receive the order information you need to achieve your business objectives. Our convenient J&J Customer Connect portal lets you:



- Enjoy 24/7 ordering & account access
- Check order status & track shipments
- View purchase history & process returns
- View your order, even if placed through EDI, phone or fax
- Easily access invoices & delivery notes
- Access backorder status by running reports or receiving daily/weekly email updates
- View product availability dates
- Contact customer support...and more!

Contents

04 Access

Endopath XCEL Trocars
Endopath Insufflation Needles
Endopath Dextrus HALS Products

10 Endo Devices

Endopath Endo Instruments
Lapra-Ty Absorbable Suture Products
Endopath Specimen Retrieval Bag
Endoloop Ligatures
Laparoscopic Needle Holder

13 Energy

Harmonic Ultrasonic Shears
Harmonic Ultrasonic Blade
Enseal Advance Bipolar Tissue Sealers
Capital Equipment & Accessories
Powerstar Bipolar Scissors
Endopath Probe Plus

17 Ligation

Endoscopic Multiple Clip Appliers
Open Multiple Clip Appliers
Endoscopic Single Clip Appliers and Clips
Open Single Clip Appliers and Clips
Absolok Absorbable Clip Appliers and Clips

17 Stapling

Echelon Powered Endocutters
Linear Cutters
Staple Line Reinforcement
Linear Cutters
Powered Circular Staplers
TX Linear Staplers
Echelon Contour Curved Cutter
Skin Staplers

17 Index

Access

Access

BLADELESS

5mm ENDOPATH XCEL Bladeless Trocars with Optiview™ Technology



Product Code	Sleeve	Shaft Length	Quantity
2B5ST	Stability	75 mm	6
2B5LT	Stability	100 mm	6
2B5XT	Stability	150 mm	6

5mm ENDOPATH XCEL Bladeless Trocars



Product Code	Sleeve	Shaft Length	Quantity
B5ST	Stability	75 mm	6
B5LT	Stability	100 mm	6
B5XT	Stability	150 mm	6

8mm ENDOPATH XCEL Bladeless Trocar



Product Code	Sleeve	Shaft Length	Quantity
B8LT	Stability	100 mm	6

11mm ENDOPATH XCEL Bladeless Trocars



Product Code	Sleeve	Shaft Length	Quantity
B11LT	Stability	100 mm	6
B11LP	Smooth	100 mm	6

11mm ENDOPATH XCEL Bladeless Trocars with Handle



Product Code	Sleeve	Shaft Length	Quantity
B11LTH	Stability	100 mm	6
B11LPH	Smooth	100 mm	6

12mm ENDOPATH XCEL Bladeless Trocars



Product Code	Sleeve	Shaft Length	Quantity
B12SRT	Stability	75 mm	6
B12LT	Stability	100 mm	6
B12XT	Stability	150 mm	6
B12LP	Smooth	100 mm	6

12mm ENDOPATH XCEL Bladeless Trocars with OPTIVIEW™ Technology



Product Code	Sleeve	Shaft Length	Quantity
2B12LT	Stability	100 mm	6
2B12XT	Stability	150 mm	6

12mm ENDOPATH XCEL Bladeless Trocars with Handle



Product Code	Sleeve	Shaft Length	Quantity
B12LTH	Stability	100 mm	6
B12LPH	Smooth	100 mm	6

15mm ENDOPATH XCEL Bladeless Trocar



Product Code	Sleeve	Shaft Length	Quantity
B15LT	Stability	100 mm	6

DILATING TIP

5mm Endopath XCEL Dilating Tip Trocars with Optiview™ Technology




Product Code	Sleeve	Shaft Length	Quantity
2D5ST	Stability	75 mm	6
2D5LT	Stability	100 mm	6

5mm ENDOPATH XCEL Dilating Tip Trocars




Product Code	Sleeve	Shaft Length	Quantity
D5ST	Stability	75 mm	6
D5LT	Stability	100 mm	6

11mm ENDOPATH XCEL Dilating Tip Trocars



Product Code	Sleeve	Shaft Length	Quantity
D11LT	Stability	100 mm	6

12mm ENDOPATH XCEL Dilating Tip Trocars with OPTIVIEW™ Technology



Product Code	Sleeve	Shaft Length	Quantity
2D12LT	Stability	100 mm	6
2D12XT	Stability	150 mm	6

12mm ENDOPATH XCEL Dilating Tip Trocars



Product Code	Sleeve	Shaft Length	Quantity
D12LT	Stability	100 mm	6
D12XT	Stability	150 mm	6

SLEEVES

5mm ENDOPATH XCEL Universal Sleeve with OPTIVIEW™ Technology

Can be used with Both XCEL™ Bladeless and Dilating Tip Trocar



Product Code	Sleeve	Shaft Length	Quantity
2CB5ST	Stability	75 mm	6
2CB5LT	Stability	100 mm	6

5mm ENDOPATH XCEL Universal Sleeve

Can be used with Both XCEL™ Bladeless and Dilating Tip Trocar



Product Code	Sleeve	Shaft Length	Quantity
CB5LT	Stability	100 mm	6
CB5ST	Stability	75 mm	6

11mm ENDOPATH XCEL Universal Sleeve



Product Code	Sleeve	Shaft Length	Quantity
CB11LT	Stability	100 mm	6

12mm ENDOPATH XCEL Universal Sleeve with OPTIVIEW™ Technology



Product Code	Sleeve	Shaft Length	Quantity
2CB12LT	Stability	100 mm	6

12mm ENDOPATH XCEL Universal Sleeve



Product Code	Sleeve	Shaft Length	Quantity
CB12LT	Stability	100 mm	6

BLUNT TIP

12mm ENDOPATH XCEL Blunt Tip Trocar with OPTIVIEW™ Technology



Product Code	Shaft Length	Quantity
2H12LP	100 mm	6

12mm ENDOPATH XCEL Blunt Tip Trocar



Product Code	Shaft Length	Quantity
H12LP	100 mm	6

NEEDLES

ENDOPATH Insufflation Needles Pneumoperitoneum



Product Code	Shaft Length	Quantity
PN120	120 mm	12
PN150	150 mm	12

ENDOPATH Insufflation Needle Ultra Veress



Product Code	Sleeve	Shaft Length	Quantity
UV120	Stability	100 mm	12

HAND-ASSISTED LAPAROSCOPY (HALS)

ENDOPATH DEXTRUS HALS Retractors



Product Code	Description	Shaft Length	Quantity
FLR01	Fixed Length Access Retractor - Small	<4 cm	6
FLR02	Fixed Length Access Retractor - Medium	4-7 cm	6
FLR03	Fixed Length Access Retractor - Large	>7 cm	6

ENDOPATH DEXTRUS HALS Products Retractor sold separately



Product Code	Description	Quantity
HAP02	Seal Cap Assembly with Accessories	6

Endo Devices

Endo Devices

ENDO INSTRUMENTS

ENDOPATH 5mm Instruments



Product Code	Shaft Length	Quantity
5BB	Babcock, 5mm	6
5DCD	Dissector, Curved, with Monopolar Cautery, 5mm	6
5DCS	Scissors, Curved, with Monopolar Cautery, 5mm	6
5DSG	Grasper, 5mm	6

ENDOPATH 5mm Instruments – Blunt Tip




Product Code	Description	Quantity
BTD05	Dissector, Blunt Tip, 5mm, 3 per pouch, 12 Pouches per unit	12

ENDOPATH 10mm Instruments



Product Code	Description	Quantity
10AG	Anvil Grasper 10 mm	6
10BB	Babcock 10 mm	6

ENDOPATH 10mm Instruments – Blunt Tip



Product Code	Description	Quantity
BCD10	Dissector, Blunt Cherry, 10mm, 3 per Pouch, 12 Pouches per Unit	12

ENDOSUTURE



LAPRA-TY Absorbable Suture Products

Product Code	Description	Quantity
MIC4030	LAPRA-TY® Absorbable Suture Clip Applier, Reusable, for use with 11mm or 12mm Diameter Trocar	1



LAPRA-TY Absorbable Suture Products

Product Code	Description	Quantity
MIC4032	LAPRA-TY® Absorbable Suture Clips for use with 2-0, 3-0, 4-0 Coated VICRYL Suture, 2 Clips per Cartridge	6
MIC4034	LAPRA-TY® Absorbable Suture Clips for use with 2-0, 3-0, 4-0 Coated VICRYL Suture, 4 Clips per Cartridge	6
MIC4036	LAPRA-TY® Absorbable Suture Clips for use with 2-0, 3-0, 4-0 Coated VICRYL Suture, 6 Clips per Cartridge	6


POUCH

ENDOPATH Specimen Retrieval Bags



Product Code	Description	Shaft	Quantity
POUCH	4"x6" 224ml, 10 mm	26 cm	6

ENDOLOOP Ligature



Product Code	Description	Quantity
EJ10C	ENDOLOOP® Coated VICRYL Ligature, 0, 50 cm	12
EJ11	ENDOLOOP®, VICRYL Loop Ligature, 0, 50 cm, non-needled	12
EZ11	ENDOLOOP®, PDS Loop Ligature, 0, 50 cm, non-needled	12
MIC108G	ENDOLOOP, PDS II Loop Ligature, 3-0	12
MIC109G	ENDOLOOP, PDS II Loop Ligature, 2-0	12
MIC110G	ENDOLOOP, PDS II Loop Ligature, 0	12
MIC111G	ENDOLOOP, PDS II, Loop Ligature, Large, 1, 70 cm	12
MIC113G	ENDOLOOP, PDS II Loop Ligature, Large, 1	12

Laparoscopic Needle Holder

Product Code	Description	Quantity
E705R	ENDOPATH, Needle Holder, Reusable, 5 mm, 30 cm Length	1
SRNH1	Self Righting Needle Holder, Reusable, 5 mm, 30 cm Length	1

Energy

Energy

All ENSEAL and HARMONIC devices are compatible with the Ethicon Endo-Surgery Generator G11 (code: GEN11).

SHEARS



HARMONIC 1100 Shears with Adaptive Tissue Technology

Product Code	Tip	Shaft Diameter	Shaft Length	Handpiece	Quantity
HAR1120	Curved Shears	5 mm	20 cm	Integrated	6
HAR1136	Curved Shears	5 mm	36 cm	Integrated	6



HARMONIC 700 Shears with Adaptive Tissue Technology

Product Code	Tip	Shaft Diameter	Shaft Length	Handpiece	Quantity
HAR723	Curved Shears	5 mm	23 cm	HARHPGR	6
HAR736	Curved Shears	5 mm	36 cm	HARHPGR	6
HAR745	Curved Shears	5 mm	45 cm	HARHPGR	6



HARMONIC ACE + Shears with Adaptive Tissue Technology

Product Code	Tip	Shaft Diameter	Shaft Length	Handpiece	Quantity
HAR23	Curved Shears	5 mm	23 cm	HARHPGR	6
HAR36	Curved Shears	5 mm	36 cm	HARHPGR	6



HARMONIC FOCUS + Shears with Adaptive Tissue Technology

Product Code	Tip	Shaft Length	Handpiece	Quantity
HAR9F	Curved Shears	9 mm	HARHPBL	6
HAR17F	Curved Shears	17 mm	HARHPBL	6

BLADES



HARMONIC Laparoscopic Blades

Product Code	Tip	Shaft Length	Handpiece	Quantity
HDH05	Hook	5 mm	HARHPGR	6

TISSUE SEALERS

All ENSEAL and HARMONIC devices are compatible with the Ethicon Endo-Surgery Generator G11 (code: GEN 11).



Product Code	Shaft Diameter	Shaft Length	Quantity
NSLX125C	5 mm	25 cm	3
NSLX137C	5 mm	37 cm	3
NSLX145C	5 mm	45 cm	3



Product Code	Shaft Diameter	Shaft Length	Quantity
NSLX125S	5 mm	25 cm	3
NSLX137S	5 mm	37 cm	3



Product Code	Shaft Diameter	Shaft Length	Quantity
NSLX120L	13 mm	20 cm	6

ENSEAL G2 Articulating Curved



Product Code	Shaft Diameter	Shaft Length	Quantity
NSLG2C35A	5 mm	35 cm	6
NSLG2C45A	5 mm	45 cm	6

ENSEAL G2 Articulating Straight



Product Code	Shaft Diameter	Shaft Length	Quantity
NSLG2S35A	5 mm	35 cm	6
NSLG2S45A	5 mm	45 cm	6

CAPITAL

EES Generator



Product Code	Description	Quantity
GEN11	EES Generator compatible with all HARMONIC™ and ENSEAL™ DEVICES	1

Generator Cart



Product Code	Description	Quantity
CRT11	Generator Cart	1

Foot Pedal



Product Code	Description	Quantity
FSW11	Generator Foot Switch	1

HANDPIECES AND ACCESSORIES



Product Code	Description	Quantity
HARHPGR	HARMONIC™ Grey Hand Piece	1
HARHPBL	HARMONIC™ Blue Hand Piece	1

Accessories













Product Code	Shaft Diameter	Quantity
HST01	Sterilization Tray	1
HST02	Test Tip for HARHPGR	1
TLB01	Torque Lock Blade Wrench	2
TTBLUE	Test Tip for HARHPBL	1
HSA08	Disposable Hand Switch Adaptor	3
GEN11VK	GEN11 Verification Key	1
HP54CTA	HARMONIC™ Hand Piece Clean Tool	1

OTHER ENERGY

POWERSTAR Bipolar Scissors

Product Code **Description**

BP100	PowerStar Tenotomy 110 mm (4 1/2') curved	
BP120	PowerStar Tenotomy 175 mm curved	
BP310	PowerStar Metzenbaum Fine 140 mm curved	
BP320	PowerStar Metzenbaum Fine 175 mm (7') curved	
BP340	PowerStar Metzenbaum Fine 230 mm (9') curved	
BP360	PowerStar Metzenbaum Fine 280 mm (11') curved	
BP520	PowerStar Metzenbaum 175 mm (7') curved	
BP540	PowerStar Metzenbaum 230 mm (9') curved	
BP560	PowerStar Metzenbaum 280 mm (11') curved	
BP720	PowerStar Mayo 170 mm (6 3/4') curved	



ENDOPATH PROBE PLUS II Handles

Product Code	Control	Shaft Length	Quantity
EPH02	Hand	Pistol	6
EPH04	Hand	Pencil	6



ENDOPATH PROBE PLUS II Hollow Tip Electrodes

Product Code	Electrode	Shaft Diameter	Shaft Diameter	Quantity
EPS01	Hook	5 cm	34 cm	6
EPS02	Spatula	5 cm	34 cm	6
EPS03	Right Angle	5 cm	34 cm	6
EPS04	Curved Disector	5 cm	34 cm	6
EPS06	Spatula	5 cm	29 cm	6
EPS07	Right Angle	5 cm	29 cm	6



Endopath Active Cords

Product Code	Description	Pins	Quantity
EBC02	Reusable Coaxial Pin Connector Bipolar Cord	1	1
EBC05	Active Return Cord Fixed Pin Connector	2	1

Ligation

Ligation

ENDOSCOPIC MULTIPLE CLIP APPLIERS

Ligamax™ 5 Rotating Multiple Clip Applier



Product Code	Diameter	Shaft Length	Clip Size	No. of Clips	Quantity
EL5ML	5 mm	33 cm	Medium/Large	15	3

Ligacip™ Rotating Multiple Clip Applier



Product Code	Diameter	Shaft Length	Clip Size	No. of Clips	Quantity
ER320	10 mm	28,9 cm	Medium/Large	20	3



Product Code	Diameter	Shaft Length	Clip Size	No. of Clips	Quantity
ER420	12 mm	34,1 cm	Large	20	3

OPEN MULTIPLE CLIP APPLIERS

Ligacip™ Multiple Clip Appliers



Product Code	Length	Clip Size	No. of Clips	Quantity
MCS20	23.8 cm	Small	20	6
MSM20	23.8 cm	Medium	20	6
MCM20	29.2 cm	Medium	20	6
MCM30	29.2 cm	Medium	30	6
MCL20	33.7 cm	Large	20	6

ENDOSCOPIC SINGLE CLIP APPLIERS

Ligacip™ Multi-Patient Use Single Clip Appliers



Product Code	Clip Size	Quantity
LC2010Z	Medium	1
LC3010Z	Medium/Large	1
LC4010Z	Large	1

OPEN SINGLE CLIP APPLIERS

Ligacip™ Multi-Patient Use Single Clip Appliers – Blue



Product Code	Handle Color	Length	Clip Size	Quantity
LX105	Blue	14.6 cm	Small	1
LX107	Blue	19.1 cm	Small	1
LX110	Blue	26.7 cm	Small	1

Ligacip™ Multi-Patient Use Single Clip Appliers – Silver



Product Code	Description	Handle Color	Length	Clip Size	Quantity
LX205	-----	Silver	14.6 cm	Medium	1
LX207	-----	Silver	19.1 cm	Medium	1
LX210	-----	Silver	26.7 cm	Medium	1
LX220	Angled Jaw	Silver	26.7 cm	Medium	1

Ligaclip™ Multi-Patient Use Single Clip Appliers – Green



Product Code	Description	Handle Color	Length	Clip Size	Quantity
LC307	-----	Green	19.1 cm	Medium/Large	1
LC310	-----	Green	26.7 cm	Medium/Large	1
LC320	Angled Jaw	Green	26.7 cm	Medium/Large	1

Ligaclip™ Multi-Patient Use Single Clip Appliers – Yellow



Product Code	Description	Handle Color	Length	Clip Size	Quantity
LC407	-----	Yellow	19.1 cm	Large	1
LC410	-----	Yellow	26.7 cm	Large	1
LC420	Angled Jaw	Yellow	26.7 cm	Large	1

Ligaclip™ Ligating Clip Cartridge Base



Product Code	Combatible With	Quantity
LC800	LIGACLIP™ Ligating Clip Cartridge	1

SINGLE CLIPS



Ligacip™ EXTRA Titanium Ligating Clips

For use with both Endoscopic and Open/Conventional Single Clip Reusable Appliers

Product Code	Color	Clip Size	No. of Clips	Quantity
LT100	Blue	Small	6	36
LT102	Blue	Small	20	15
LT200	Silver	Medium	6	36
LT202	Silver	Medium	20	15
LT300	Green	Medium/Large	6	18
LT400	Yellow	Large	6	18

ABSOLOK ABSORBABLE CLIP APPLIERS AND CLIPS



ABSOLOK Endoscopic Reusable Single Ligating Clip Applier

Product Code	Diameter	Clip Size	Color Code	Quantity
AP2010S	10 mm	Medium	● Silver	1
AP3010S	10 mm	Medium/Large	● Green	1
AP4010S	10 mm	Large	● Yellow	1



ABSOLOK Open Reusable Single Ligating Clip Applier

Product Code	Length	Clip Size	Color Code	Quantity
AC105	14,5 cm	Small	● Blue	1
AC107	18,5 cm	Small	● Blue	1
AC207	18,5 cm	Medium	● Silver	1
AC210	27,5 cm	Medium	● Silver	1
AC307	18,5 cm	Medium/Large	● Green	1
AC310	27,5 cm	Medium/Large	● Green	1
AC410	27,5 cm	Large	● Yellow	1



ABSOLOK Absorbable Single Clips

Product Code	Color Code	Clip Size	No. of Clips	Quantity
AP100	● Blue	Small	10	6
AP200	● Silver	Medium	10	6
AP201	● Silver	Medium	5	6
AP300	● Green	Medium/Large	10	6
AP301	● Green	Medium/Large	5	6
AP306	● Green	Medium/Large	6	6
AP400	● Yellow	Large	10	6
AP401	● Yellow	Large	5	6

Stapling

Stapling

SURGICAL STAPLING

Powered ECHELON 3000 60mm Endocutters

Reloads Sold Separately



Product Code	Articulation	Powered	Shaft Length	Quantity
ECH60C	55	Yes	280 mm	3
ECH60S	55	Yes	340 mm	3
ECH60L	55	Yes	440 mm	3

Powered Plus ECHELON FLEX 60mm Endocutters

Reloads Sold Separately



Product Code	Articulation	Powered	Shaft Length	Quantity
PCEE60A	45	Yes	280 mm	3
PSEE60A	45	Yes	340 mm	3
PLEE60A	45	Yes	440 mm	3

ECHELON FLEX 60mm Endocutters

Reloads Sold Separately



Product Code	Articulation	Powered	Shaft Length	Quantity
SC60A	45	No	280 mm	3
EC60A	45	No	340 mm	3
LONG60A	45	No	440 mm	3

Powered ECHELON 3000 45mm Endocutters

Reloads Sold Separately



Product Code	Articulation	Powered	Shaft Length	Quantity
ECH45C	55	Yes	280	3
ECH45S	55	Yes	340	3
ECH45L	55	Yes	440	3

Powered Plus ECHELON FLEX 45mm Endocutters

Reloads Sold Separately



Product Code	Articulation	Powered	Shaft Length	Quantity
PCEE45A	45	Yes	280	3
PSEE45A	45	Yes	340	3
PLEE45A	45	Yes	440	3

ECHELON FLEX 45mm Endocutters

Reloads Sold Separately



Product Code	Articulation	Powered	Shaft Length	Quantity
SC45A	45	No	280	3
EC45A	45	No	340	3
EC45AL	45	No	440	3



ECHELON 60mm Reloads with GST

Product Code	Color	Rows	Tissue Type	Open Staple Height	Closed Staple Height	Quantity
GST60W	White	6	Vascular	2.6 mm	1.0 mm	12
GST60B	Blue	6	Standard	3.6 mm	1.5 mm	12
GST60D	Gold	6	Standard/Thick	3.8 mm	1.8 mm	12
GST60G	Green	6	Thick	4.1 mm	2.0 mm	12
GST60T	Black	6	Extra Thick	4.2 mm	2.3 mm	12

ECHELON 45mm Reloads with GST

Product Code	Color	Rows	Tissue Type	Open Staple Height	Closed Staple Height	Quantity
GST45W	White	6	Vascular	2.6 mm	1.0 mm	12
GST45B	Blue	6	Standard	3.6 mm	1.5 mm	12
GST45D	Gold	6	Standard/Thick	3.8 mm	1.8 mm	12
GST45G	Green	6	Thick	4.1 mm	2.0 mm	12
GST45T	Black	6	Extra Thick	4.2 mm	2.3 mm	12

Staple Line Reinforcement



Product Code	Compatible With	Quantity
ECH60R	ECHELON 60mm GST Reloads	6

ECHELON FLEX Powered Vascular Stapler

Reloads Sold Separately



Product Code	Articulation	Powered	Shaft Length	Quantity
PVE35A	Yes	Yes	320 mm	3

ENDOPATH ECHELON Vascular White Reload 35mm for Advanced Placement Tip



Product Code	Color	Rows	Tissue Type	Open Staple Height	Closed Staple Height	Quantity
VASECR35	<input type="radio"/> White	4	Vascular	2.5 mm	1.0 mm	12

ETS Flex 45mm Endocutters

Reloads Sold Separately



Product Code	Articulation	Powered	Shaft Length	Quantity
ATS45	Yes	No	340 mm	3

ETS 45mm Reloads

For use with all ETS 45mm Endoscopic Linear Cutters ATS45



Product Code	Color	Rows	Tissue Type	Open Staple Height	Closed Staple Height	Quantity
TR45W	<input type="radio"/> White	6	Vascular	2.5 mm	1.0 mm	12
6R45B	<input checked="" type="radio"/> Blue	6	Standard	3.5 mm	1.5 mm	12

LINEAR CUTTERS

ECHELON LINEAR CUTTER Stapler 60mm

Reloads Sold Separately



Product Code	Description	Jaw Length	Quantity
GLC60	ECHELON LINEAR Cutter 60mm	60 mm	3

ECHELON LINEAR CUTTER Reload

Reloads Sold Separately



Product Code	Color	Rows	Tissue Range	Open Staple Height	Closed Staple Height	Quantity
GLCR60B	Blue	4	1,5mm - 2,4mm	3.7 mm	1.5 mm	12
GLCR60G	Green	4	2,0mm - 3,3mm	4.5 mm	2.0 mm	12

ECHELON LINEAR CUTTER Stapler 80mm

Reloads Sold Separately



Product Code	Description	Jaw Length	Quantity
GLC80	ECHELON LINEAR Cutter 80mm	80 mm	3

ECHELON LINEAR CUTTER Reload

Reloads Sold Separately



Product Code	Color	Rows	Tissue Range	Open Staple Height	Closed Staple Height	Quantity
GLCR80B	Blue	4	1,5mm - 2,4mm	3.7 mm	1.5 mm	12
GLCR80G	Green	4	2,0mm - 3,3mm	4.5 mm	2.0 mm	12

ECHELON LINEAR CUTTER Stapler 100mm

Reloads Sold Separately





Product Code	Description	Jaw Length	Quantity
GLC100	ECHELON LINEAR Cutter 100mm	100 mm	3

ECHELON LINEAR CUTTER Reload

Reloads Sold Separately



Product Code	Color	Rows	Tissue Range	Open Staple Height	Closed Staple Height	Quantity
GLCR100B	 Blue	4	1,5 mm - 2,4 mm	3.7 mm	1.5 mm	12
GLCR100G	 Green	4	2,0 mm - 3,3 mm	4.5 mm	2.0 mm	12

Ethicon Linear Cutter 55, Selectable Staple Height




Reloads Sold Separately



Product Code	Selectable Colour	Tissue Type	Quantity
NTLC55	  	Selectable	3

Ethicon Linear Cutter 55, Selectable Staple Height – Reload



Product Code	Color	Rows	Tissue Range	Open Staple Height	Closed Staple Height	Quantity
SR55	 Blue	6	Standard	4.3 mm	1.5 mm	12
	 Gold	6	Standard/Thick	4.3 mm	1.8 mm	12
	 Green	6	Thick	4.3 mm	2.0 mm	12

Ethicon Linear Cutter 75, Selectable Staple Height

Reloads Sold Separately



Product Code	Selectable Colour	Tissue Type	Quantity
NTLC75	<div><div></div><div></div><div></div></div>	Selectable	3

Ethicon Linear Cutter 75, Selectable Staple Height – Reload



Product Code	Color	Rows	Tissue Range	Open Staple Height	Closed Staple Height	Quantity
SR75	<div></div> Blue	6	Standard	4.3 mm	1.5 mm	12
	<div></div> Gold	6	Standard/Thick	4.3 mm	1.8 mm	12
	<div></div> Green	6	Thick	4.3 mm	2.0 mm	12

55mm Linear Cutters

PROXIMATE 55mm Linear Cutters with Safety Lockout



Product Code	Color	Tissue Type	Reload Included	Open Staple Height	Closed Staple Height	Quantity
TVC55	<div></div> White	Vascular	TVR55	2.5 mm	1.0 mm	3
TLC55	<div></div> Blue	Standard	TCR55	3.85 mm	1.5 mm	3
TCT55	<div></div> Green	Thick	TRT55	4.5 mm	2.0 mm	3

55mm Linear Cutter Reloads

For use with all 55mm Linear Cutters (TVC55, TLC55, TCT55)






Product Code	Color	Rows	Tissue Range	Open Staple Height	Closed Staple Height	Quantity
TVR55	<div></div> White	4	Vascular	2.5 mm	1.0 mm	12
TCR55	<div></div> Blue	4	Standard	3.85 mm	1.5 mm	12
TRT55	<div></div> Green	4	Thick	4.50 mm	2.0 mm	12

75mm Linear Cutters

PROXIMATE 55mm Linear Cutters with Safety Lockout






Product Code	Color	Tissue Type	Reload Included	Open Staple Height	Closed Staple Height	Quantity
TLC75	 Blue	Standard	TCR75	3.85 mm	1.0 mm	3
TCD75	 Gold	Standard/Thick	TRD75	4.2 mm	1.8 mm	3
TCT75	 Green	Thick	TRT75	4.5 mm	2.0 mm	3

75mm Linear Cutter Reloads

For use with all 75mm Linear Cutters (TLC75, TCD75, TCT75)





Product Code	Color	Rows	Tissue Type	Open Staple Height	Closed Staple Height	Quantity
TCR75	 Blue	4	Standard	3.85 mm	1.0 mm	12
TRD75	 Gold	4	Standard/Thick	4.2 mm	1.5 mm	12
TRT75	 Green	4	Thick	4.5 mm	2.0 mm	12

100mm Linear Cutters

PROXIMATE® 100mm Linear Cutters with Safety Lockout





Product Code	Color	Tissue Type	Reload Included	Open Staple Height	Closed Staple Height	Quantity
TLC10	 Blue	Standard	TCR10	3.85 mm	1.5 mm	3
TCT10	 Green	Thick	TRT10	4.5 mm	2.0 mm	3

100mm Linear Cutter Reloads

For use with all 100mm Linear Cutters (TLC10, TCT10)



Product Code	Color	Rows	Tissue Type	Open Staple Height	Closed Staple Height	Quantity
TCR10	 Blue	4	Standard	3.85 mm	1.5 mm	12
TRT10	 Green	4	Thick	4.5 mm	2.0 mm	12

CIRCULAR STAPLERS

Endoscopic Curved Intraluminal Staplers



Product Code	Lumen	Shaft Length	Head Diameter	Open Staple Height	Closed Staple Height	Quantity
ECS21B	12.4 mm	XL	21 mm	5.2 mm	1.5mm - 2.2mm Adjustable	3
ECS25B	16.5 mm	XL	25 mm	5.2 mm	1.5mm - 2.2mm Adjustable	3
ECS29B	20.4 mm	XL	29 mm	5.2 mm	1.5mm - 2.2mm Adjustable	3
ECS33B	24.4 mm	XL	33 mm	5.2 mm	1.5mm - 2.2mm Adjustable	3

Curved Intraluminal Staplers



Product Code	Lumen	Shaft Length	Head Diameter	Open Staple Height	Closed Staple Height	Quantity
CDH21B	12.4 mm	Standard	21 mm	5.2 mm	1.5mm - 2.2mm Adjustable	3
CDH25B	16.5 mm	Standard	25 mm	5.2 mm	1.5mm - 2.2mm Adjustable	3
CDH29B	20.4 mm	Standard	29 mm	5.2 mm	1.5mm - 2.2mm Adjustable	3
CDH33B	24.4 mm	Standard	33 mm	5.2 mm	1.5mm - 2.2mm Adjustable	3

Echelon Circular Powered Staplers



Product Code	Lumen	Shaft Length	Head Diameter	Open Staple Height	Closed Staple Height	Quantity
CDH23P	14.6 mm	24 cm	23 mm	5.2 mm	1.5mm - 2.2mm	3
CDH25P	16.5 mm	24 cm	25 mm	5.2 mm	1.5mm - 2.2mm	3
CDH29P	20.4 mm	24 cm	29 mm	5.2 mm	1.5mm - 2.2mm	3
CDH31P	22.5 mm	24 cm	31 mm	5.2 mm	1.5mm - 2.2mm	3

Circular Sizers



Product Code	Description	Quantity
CS23	Bowel sizers for the ENDOPATH ILS Curved Staplers, 23 mm	1
CS31	Bowel sizers for the ENDOPATH ILS Curved Staplers, 31 mm	1
CSS	Bowel sizers for the ENDOPATH ILS Curved Staplers, the set includes one each of a 25 mm, 29 mm, and 33 mm sizer	1

Reusable Purse String Clamp



Product Code	Description	Quantity
EH40	----	1

PROCEDURE STAPLERS

PPH Procedure Set



Product Code	Description	Quantity
PPH03	Includes Hemorrhoidal Circular Stapler 33 mm, suture threader, circular anal dilator, and purse string anoscope	1

LINEAR STAPLERS

TX 30mm Linear Stapler

PROXIMATE® TX 30mm Reloadable Linear Stapler and Reload



Product Code	Color	Tissue Type	Reload Included	Open Staple Height	Closed Staple Height	Quantity
TX30V	White	Vascular	XR30V	2.5 mm	1.0 mm	3
XR30V	White	Vascular	—	2.5 mm	1.0 mm	12

TX 30mm Linear Stapler


PROXIMATE® TX 30mm Reloadable Linear Stapler



Product Code	Color	Tissue Type	Reload Included	Open Staple Height	Closed Staple Height	Quantity
TX30B	Blue	Standard	XR30B	3.5 mm	1.0 mm	3
TX30G	Green	Thick	XR30G	4.8 mm	1.0 mm	3

TX 30mm Linear Stapler Reloads

For use with (TX30B, TX30G)



Product Code	Color	Tissue Type	Open Staple Height	Closed Staple Height	Quantity
XR30B	Blue	Standard	3.5 mm	1.0 mm	12
XR30G	Green	Thick	4.8 mm	1.0 mm	12

TX 60mm Linear Stapler

PROXIMATE® TX 60mm Reloadable Linear Stapler



Product Code	Color	Tissue Type	Reload Included	Open Staple Height	Closed Staple Height	Quantity
TX60B	Blue	Standard	XR60B	3.5 mm	1.5 mm	3
TR60G	Green	Thick	XR60G	4.8 mm	2.0 mm	3

TX 30mm Linear Stapler Reloads

For use with (TX60B, TX60G)



Product Code	Color	Tissue Type	Open Staple Height	Closed Staple Height	Quantity
XR60B	Blue	Standard	3.5 mm	1.5 mm	12
XR60G	Green	Thick	4.8 mm	2.0 mm	12

CURVED CUTTER STAPLERS

ECHELON Contour 40mm Linear Staplers



Product Code	Color	Reload Included	Open Staple Height	Closed Staple Height	Quantity
GCS40B	Blue	GCR40B	3.5 mm	1.5 mm	3
GCS40G	Green	GCR40G	4.7 mm	2.0 mm	3

ECHELON Contour 40mm Linear Stapler Reloads

For use with (GCS40B, GCS40G)



Product Code	Color	Tissue Type	Open Staple Height	Closed Staple Height	Quantity
GCR40B	Blue	Standard	3.5 mm	1.5 mm	6
GCR40G	Green	Thick	4.7 mm	2.0 mm	6

SKIN STAPLERS

Rotating Head Skin Staplers

PROXIMATE™ Rotating Head Skin Staplers with Stainless Steel Staples



Product Code	Size	Tissue Type	Quantity
PRR35	Regular	35	6
PRW35	Wide	35	6

Fixed Head Skin Staplers

PROXIMATE™ Plus MD Fixed Head Skin Staplers with Multi-Directional Release



Product Code	Size	Tissue Type	Quantity
PMW35	Wide	35	6

Fixed Head Skin Staplers Pistol Grip

PROXIMATE™ Plus MD Fixed Head Skin Staplers with Multi-Directional Release



Product Code	Size	Tissue Type	Quantity
PXW35	Wide	35	6

Staple Extractor



Product Code	Size	Tissue Type	Quantity
PSX	Squeeze Handle Staple Extractor	35	12

Index

Index

CODE	PAGE
10AG	16
10BB	16
2B12LT	9
2B12XT	9
2B5LT	7
2B5ST	7
2B5XT	7
2CB12LT	11
2CB5LT	11
2CB5ST	11
2D12LT	10
2D12XT	10
2D5LT	9
2D5ST	9
2H12LP	12
5BB	16
5DCD	16
5DCS	16
5DSG	16
6R45B	39
AC105	34
AC107	34
AC207	34
AC210	34
AC307	34
AC310	34
AC410	34
AP100	34
AP200	34
AP201	34
AP2010S	34
AP300	34
AP301	34
AP3010S	34
AP306	34
AP400	34
AP401	34
AP4010S	34
ATS45	39

CODE	PAGE
B11LP	8
B11LPH	8
B11LT	8
B11LTH	8
B12LP	8
B12LPH	9
B12LT	8
B12LTH	9
B12SRT	8
B12XT	8
B15LT	9
B5LT	7
B5ST	7
B5XT	7
B8LT	7
BCD10	17
BP100	26
BP120	26
BP310	26
BP320	26
BP340	26
BP360	26
BP520	26
BP540	26
BP560	26
BP720	26
BTD05	16
CB11LT	11
CB12LT	12
CB5LT	11
CB5ST	11
CDH21B	44
CDH23P	44
CDH25B	44
CDH25P	44
CDH29B	44
CDH29P	44
CDH31P	44
CDH33B	44

CODE	PAGE
CRT11	24
CS23	45
CS31	45
CSS	45
D11LT	10
D12LT	10
D12XT	10
D5LT	10
D5ST	10
E705R	18
EBC02	28
EBC05	28
EC45A	37
EC45AL	37
EC60A	36
ECH45C	37
ECH45L	37
ECH45S	37
ECH60C	36
ECH60L	36
ECH60R	38
ECH60S	36
ECS21B	44
ECS25B	44
ECS29B	44
ECS33B	44
EH40	45
EJ10C	18
EJ11	18
EL5ML	30
EPH02	27
EPH04	27
EPS01	27
EPS02	27
EPS03	27
EPS04	27
EPS06	27
EPS07	27
EPS10	27

CODE	PAGE
ER320	30
ER420	30
EZ11	18
FLR01	14
FLR02	14
FLR03	14
FSW11	24
GCS40B	47
GCS40G	47
GCR40B	47
GCR40G	47
GEN11	24
GEN11VK	25
GLC60	40
GLC80	40
GLC100	41
GLCR60B	40
GLCR60G	40
GLCR80B	40
GLCR80G	40
GLCR100B	41
GLCR100G	41
GST45B	38
GST45D	38
GST45G	38
GST45T	38
GST45W	40
GST60B	38
GST60D	38
GST60G	38
GST60T	38
GST60W	38
H12LP	12
HAP02	14
HAR1120	20
HAR1136	20
HAR9F	21
HAR17F	21
HAR23	20

CODE	PAGE
HAR36	20
HAR723	20
HAR736	20
HAR745	20
HARHPBL	25
HARHPGR	25
HDHO5	21
HP54CTA	25
HSA08	25
HST01	25
HST02	25
LC2010Z	31
LC3010Z	31
LC4010Z	31
LC307	32
LC310	32
LC320	32
LC407	32
LC410	32
LC420	32
LC800	33
LONG60A	36
LT100	33
LT102	33
LT200	33
LT202	33
LT300	33
LT400	33
LX105	31
LX107	31
LX110	31
LX205	31
LX207	31
LX210	31
LX220	31
MCL20	30
MCM20	30
MCM30	30
MCS20	30

CODE	PAGE
MIC108G	18
MIC109G	18
MIC110G	18
MIC111G	18
MIC113G	18
MIC4030	17
MIC4032	17
MIC4034	17
MIC4036	17
MSM20	30
NSLG2C35A	22
NSLG2C45A	22
NSLG2S35A	22
NSLG2S45A	22
NSLX120L	23
NSLX125C	23
NSLX125S	22
NSLX137S	22
NSLX137C	23
NSLX145C	23
NTLC55	41
NTLC75	42
PCEE45A	37
PCEE60A	36
PLEE45A	37
PLEE60A	36
PMW35	48
PN120	13
PN150	13
POUCH	17
PPH03	45
PRR35	48
PRW35	48
PSEE45A	37
PSEE60A	36
PSX	48
PVE35A	39
PXW35	48
SC45A	37

CODE	PAGE
SC60A	36
SR55	41
SR75	42
SRNH1	18
TCD75	43
TCR10	43
TCR55	42
TCR75	43
TCT10	43
TCT55	42
TCT75	43
TLB01	25
TLC10	43
TLC55	42
TLC75	43
TR45W	39
TRD75	43
TRT10	43
TRT55	42
TRT75	42
TTBLUE	25
TVC55	42
TVR55	42
TX30B	46
TX30G	46
TX30V	46
TX60B	46
TX60G	46
UV120	13
VASECR35	39
XR30B	46
XR30G	46
XR30V	46
XR60B	46
XR60G	47

For Finland: Products are CE marked.
BSI Group The Netherlands B.V. - Notified Body nr. 2797

Please refer always to the Instructions for Use / Package Insert
that come with the device for the most current and complete
instructions.

© Johnson & Johnson Nordic AB 2025
Kolonnvägen 45 SE - 170 67 Solna, Sweden

All rights reserved.

www.jnjmedtech.com

M_SE_ETH_ALLE_394739

son & John