

ETHICON ENERGY CATALOG

2021



CONTENTS

HARMONIC™ Ultrasonic Technology	4	Smoke Evacuation Pencils	40
HARMONIC™ HD 1000i Shears	5	MEGADYNE™ Telescoping Smoke Evacuation Pencil	41
HARMONIC ACE™+7 Shears	5	Smoke Evacuation Pencils, Shrouds, and Disposables	42
HARMONIC ACE™+ Shears	5	Smoke Evacuation Accessories	43
HARMONIC FOCUS™+ Shears	5	Stainless Steel Electrodes & Pencils	44
HARMONIC FOCUS™+ Long Shears	5	Stainless Steel Blade Electrodes	45
HARMONIC SYNERGY™ Blades	6	Stainless Steel Needle Electrodes & Ball Electrodes	46
HARMONIC™ Blade	6	Disposable Pencils With Stainless Steel Electrodes	47
HARMONIC™ Hand Piece	6	Suction Irrigation Devices	48
HARMONIC™ BLUE Hand Piece	6	Endopath™ Probe Plus™ II	49
HARMONIC™ Accessories	7	Suction Coagulators	50
ENSEAL™ Advanced Bipolar Energy	8	LLETZ Loop Electrodes	52
ENSEAL™ X1 Curved Jaw	9	POWERSTAR™ Bipolar Scissors	54
ENSEAL™ X1 Large Jaw	9	ENDOPATH™ Bipolar Forceps	58
ENSEAL™ G2 Articulating Tissue Sealers	9		
ENSEAL™ G2 Curved Tissue Sealers	10		
ENSEAL™ G2 Straight Tissue Sealers	10		
ENSEAL™ Trio Tissue Sealers	11		
ENSEAL™ Round Tip Tissue Sealers	11		
HARMONIC™ & ENSEAL™ Generator	12		
ETHICON GEN11 Generator & Accessories	13		
MEGADYNE™ Core Energy	14		
Electrosurgical Generator	16		
ETHICON MEGADYNE™ Electrosurgical Generator & Accessories	17		
E-Z CLEAN™ Non-Stick Electrosurgical Electrodes	18		
Coated Blade Electrodes	19		
Coated Needle Electrodes	20		
Blunt Needle Electrodes	20		
Coated Ball Electrodes	21		
Coated Specialty Electrodes	21		
MEGA FINE™ Needle Electrodes	22		
Protective Nose Cone E-Z CLEAN™ Electrodes	23		
Disposable Laparoscopic Electrodes	24-25		
Accessories	25		
Electrosurgical Pencils with E-Z CLEAN™	26		
Disposable Pencils with E-Z CLEAN™ Electrodes	27		
Telescoping Soft Tissue Dissectors	28		
MEGADYNE™ Telescoping Soft Tissue Dissectors	29		
Reusable Patient Return Electrodes	30		
MEGA SOFT™ Reusable Return Electrodes	31		
MEGA SOFT™ Cables	32		
Disposable Patient Return Electrodes	34		
Pre-Attached Cords	35		
Reusable Cords	35		
Surgical Smoke Evacuators	36		
MEGADYNE™ Smoke Evacuator & Accessories	37		
MINI VAC™ Smoke Evacuator & Accessories	38		

THE MOST COMPREHENSIVE PORTFOLIO OF INTELLIGENT ENERGY SOLUTIONS*




HARMONIC™ ENSEAL™ MEGADYNE™

*Upon review of 2017 product catalogs and websites (Advanced Energy, Core Energy, Surgical Smoke Evacuation) for the respective manufacturers: Medtronic, Olympus, Intuitive Surgical, Conmed (121066-191128)

HARMONIC™

Ultrasonic Technology



HARMONIC™

HARMONIC™

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 HARMONIC™ HD 1000i Shears Shaft Diameter 5mm, Shaft Length 20cm, Tip Curved Shaft Diameter 5mm, Shaft Length 36cm, Tip Curved	HARHD20	6
	HARHD36	6
 HARMONIC ACE™+7 Shears Shaft Diameter 5mm, Shaft Length 23cm, Tip Curved Shaft Diameter 5mm, Shaft Length 36cm, Tip Curved Shaft Diameter 5mm, Shaft Length 45cm, Tip Curved	HARH23	6
	HARH36	6
	HARH45	6
 HARMONIC ACE™+ Shears Shaft Diameter 5mm, Shaft Length 23cm, Tip Curved Shaft Diameter 5mm, Shaft Length 36cm, Tip Curved	HAR23	6
	HAR36	6
 HARMONIC FOCUS™+ Shears Length 9cm, Tip Curved	HAR9F	6
 HARMONIC FOCUS™+ Long Shears Length 17cm, Tip Curved	HAR17F	6

HARMONIC™

HARMONIC™

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 HARMONIC SYNERGY™ Blades Shaft Length Adjustable 4-9cm, Tip Curved Shaft Length Adjustable 4-9cm, Tip Hook Shaft Length Adjustable 4-9cm, Combination Hook	SNGCB	6
	SNGHK	6
	SNGHK2	6
 HARMONIC™ Blade Shaft Diameter 5mm, Shaft Length 32cm, Tip Hook	HDH05	6
 HARMONIC™ Hand Piece Compatible with HAR23, HAR36, HARH23, HARH36, HARH45, HDH05	HP054	1
 HARMONIC™ BLUE Hand Piece Compatible with HAR9F, HAR17F, SNGCB, SNGHK, SNGHK2	HPBLUE	1


HARMONIC™

Accessories

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
Disposable hand switch adapter (only use with HP054)	HSA08	3
Sterilization tray	HST01	1
TestTip for HP054	HST02	1
Torque lock blade wrench	TLB01	2
TestTip for HPBLUE	TTBLUE	1

ENSEAL™

Advanced Bipolar Energy



ENSEAL™

ENSEAL™

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 ENSEAL™ X1 Curved Jaw Shaft Diameter 5mm, Shaft Length 25cm, Jaw Curved Shaft Diameter 5mm, Shaft Length 37cm, Jaw Curved Shaft Diameter 5mm, Shaft Length 45cm, Jaw Curved	NSLX125C	3
	NSLX137C	3
	NSLX145C	3
 ENSEAL™ X1 Large Jaw Shaft Diameter 13mm, Shaft Length 20cm, Jaw Curved	NSLX120L	6

ENSEAL™

ENSEAL™

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
		
ENSEAL™ G2 Articulating Tissue Sealers		
Shaft Diameter 5mm, Shaft Length 35cm, Jaw Curved	NSLG2C35A	6
Shaft Diameter 5mm, Shaft Length 45cm, Jaw Curved	NSLG2C45A	6
Shaft Diameter 5mm, Shaft Length 35cm, Jaw Straight	NSLG2S35A	6
Shaft Diameter 5mm, Shaft Length 45cm, Jaw Straight	NSLG2S45A	6
		
ENSEAL™ G2 Curved Tissue Sealers		
Shaft Diameter 5mm, Shaft Length 14cm, Jaw Curved	NSLG2C14	6
Shaft Diameter 5mm, Shaft Length 25cm, Jaw Curved	NSLG2C25	6
Shaft Diameter 5mm, Shaft Length 35cm, Jaw Curved	NSLG2C35	6
Shaft Diameter 5mm, Shaft Length 45cm, Jaw Curved	NSLG2C45	6
		
ENSEAL™ G2 Straight Tissue Sealers		
Shaft Diameter 5mm, Shaft Length 14cm, Jaw Straight	NSLG2S14	6
Shaft Diameter 5mm, Shaft Length 25cm, Jaw Straight	NSLG2S25	6
Shaft Diameter 5mm, Shaft Length 35cm, Jaw Straight	NSLG2S35	6
Shaft Diameter 5mm, Shaft Length 45cm, Jaw Straight	NSLG2S45	6

ENSEAL™

ENSEAL™

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 ENSEAL™ Trio Tissue Sealers Shaft Diameter 5mm, Shaft Length 14cm, Jaw Curved Shaft Diameter 5mm, Shaft Length 25cm, Jaw Curved Shaft Diameter 5mm, Shaft Length 35cm, Jaw Curved Shaft Diameter 5mm, Shaft Length 45cm, Jaw Curved	ETRIO314H	6
	ETRIO325H	6
	ETRIO335H	6
	ETRIO345H	6

 ENSEAL™ Round Tip Tissue Sealers Shaft Diameter 5mm, Shaft Length 14cm, Jaw Straight Shaft Diameter 5mm, Shaft Length 25cm, Jaw Straight Shaft Diameter 5mm, Shaft Length 35cm, Jaw Straight Shaft Diameter 5mm, Shaft Length 45cm, Jaw Straight	NSEAL514RH	6
	NSEAL525RH	6
	NSEAL535RH	6
	NSEAL545RH	6

HARMONIC™ & ENSEAL™ GENERATOR



GENERATOR

ETHICON GEN11 GENERATOR

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 ETHICON GEN11 GENERATOR	GEN11	1

Accessories

 HPO54 and HPBLUE hand pieces	HGA11	1
 ENSEAL™ Trio Tissue Sealer and ENSEAL™ RoundTip Tissue Sealer	EGA11	1
 Output verification key	GEN11VK	1
 Foot switch and cable	FSW11	1
 Generator cart	96007	1

MEGADYNE™

Core Energy



Electrosurgical Generator

E-Z CLEAN™ Non-Stick Electrosurgical Electrodes

Electrosurgical Pencils with E-Z CLEAN™

Telescoping Soft Tissue Dissectors

Reusable Patient Return Electrodes

Disposable Patient Return Electrodes

Surgical Smoke Evacuators

Smoke Evacuation Pencils

Stainless Steel Electrodes & Pencils

Suction Irrigation Devices

Suction Coagulators

LLETZ Loop Electrodes

MEGADYNE™ ELECTROSURGICAL GENERATOR



ELECTROSURGICAL GENERATOR

ETHICON MEGADYNE™ Electrosurgical Generator

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
		
ETHICON MEGADYNE™ Electrosurgical Generator	MEGEN1	1

Accessories



MEGADYNE™ Mega Cart with Top Shelf	1300SJ	1
------------------------------------	--------	---



Monopolar Footswitch (3.7m cable)	1400JJ	1
-----------------------------------	--------	---



Round Bipolar Footswitch (3.7m cable)	1459J	1
---------------------------------------	-------	---



Bipolar Footswitch (3.7m cable)	1450J	1
---------------------------------	-------	---



Monopolar Adapter	96007	1
-------------------	-------	---




Disposable Flying Lead Bipolar Cable (3.7m) 10/box	4000J	1
--	-------	---



Disposable Fixed Lead Bipolar Cable (3.7m) 10/box	4005J	1
---	-------	---

**E-Z CLEAN™
NON-STICK
ELECTROSURGICAL
ELECTRODES**



E-Z CLEAN™ NON-STICK ELECTROSURGICAL ELECTRODES

Coated Blade Electrodes

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 2.5 inch (6.35 cm), Standard Blade Electrode	0012	12
 2.5 inch (6.35 cm), Standard Blade Electrode Modified	0012M	12
 2.75 inch (7 cm), X-Long Standard Blade Electrode	0012A	12
 2.75 inch (7 cm), X-Long Standard Blade Electrode Modified	0012AM	12
 4 inch (10.2 cm), Extended Blade Electrode	0014A	12
 4 inch (10.2 cm), Extended Blade Electrode Modified	0014AM	12
 6.5 inch (16.5 cm), Extended Blade Electrode	0014	12
 6.5 inch (16.5 cm), Extended Blade Electrode Modified	0014M	12

E-Z CLEAN™ NON-STICK ELECTROSURGICAL ELECTRODES

Coated Needle Electrodes

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 2.75 inch (7 cm), Needle Electrode	0013	12
 2.75 inch (7 cm), Needle Electrode Modified	0013M	12
 4 inch (10.2 cm), Extended Needle Electrode	0016A	12
 4 inch (10.2 cm), Needle Electrode Modified	0016AM	12
 6 inch (15.2 cm), Extended Needle Electrode	0016	12
 6 inch (15.2 cm), Extended Needle Electrode Modified	0016M	12

Blunt Needle Electrodes

 2.75 inch (7 cm), Blunt Needle Electrode	0113	12
 2.75 inch (7 cm), Blunt Needle Electrode Modified	0113M	12
 4.5 inches (11.4 cm), Blunt Needle Electrode	0113A	12

E-Z CLEAN™ NON-STICK ELECTROSURGICAL ELECTRODES

Coated Ball Electrodes

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 2 inch (5.1 cm), Ball Electrode	0015	12
 5 inch (12.7 cm), Ball Electrode	0009	12

Coated Specialty Electrodes

 5.8 inch (14.8 cm), Bayonet Needle Electrode	0028	12
 5.8 inch (14.8 cm), Bayonet Needle Electrode Modified	0028M	12
 6.3 inch (16.1 cm), Bayonet Blade Electrode Modified	0029M	12

E-Z CLEAN™ NON-STICK ELECTROSURGICAL ELECTRODES

MEGA FINE™ Needle Electrodes

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 2 inch (5.1 cm), MEGA FINE™ Needle Electrode	O118	12
 2.5 inch (6.35 cm), Extended MEGA FINE™ Needle Electrode	O118A	12
 2 inch (5.1 cm), Angled MEGA FINE™ 45 Degree	O119	12
 2 inch (5.1 cm), Angled MEGA FINE™ 45 Degree 4mm	O119A	12
 2 inch (5.1 cm), Angled MEGA FINE™ 90 Degree	O120	12
 6.5 inch (16.5 cm) MEGA FINE™ Needle	O121	12

E-Z CLEAN™ NON-STICK ELECTROSURGICAL ELECTRODES

Protective Nose Cone Electrodes

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 2.5 inch (6.35 cm), Blade Electrode Modified with PTFE Insulation and Nose Cone	0012MD	24
 2.75 inch (7 cm), Blade Electrode Modified with PTFE Insulation and Nose Cone	0012AMD	24
 4 inch (10.2 cm), Blade Electrode Modified with PTFE Insulation and Nose Cone	0014AMD	24
 6.5 inch (16.5 cm), Blade Electrode Modified with PTFE Insulation and Nose Cone	0014MD	24
 2.75 inch (7 cm), Needle Electrode Modified with PTFE Insulation and Nose Cone	0013MD	24

E-Z CLEAN™ NON-STICK ELECTROSURGICAL ELECTRODES

Disposable Laparoscopic Electrodes

DESCRIPTION	STANDARD CAT# Compatible With Standard Pencils	SPLIT CAT# Compatible With Reusable Foot Control Cords	QTY PER BOX
 33 cm Standard Blade	0017		6
 33 cm Spatula	0018		6
 33 cm Curved Spatula	0018C	0018CS	6
 33 cm Curved	0019	0019S	6
 45 cm Curved	0019L	0019LS	6
 33 cm L-Hook	0020	0020S	6
 45 cm L-Hook	0020L	0020LS	6
 33 cm J-Hook	0021	0021S	6
 45 cm J-Hook	0021L	0021LS	6
 33 cm Needle	0022	0022S	6

E-Z CLEAN™ NON-STICK ELECTROSURGICAL ELECTRODES

Disposable Laparoscopic Electrodes

DESCRIPTION	STANDARD CAT# Compatible With Standard Pencils	SPLIT CAT# Compatible With Reusable Foot Control Cords	QTY PER BOX
 33 cm Ball	0024		6
 33 cm L-Wire	0100	0100S	6
 45 cm L-Wire	0100L	0100LS	6

Accessory

DESCRIPTION	STANDARD CAT#	QTY PER BOX
 Reusable Foot Control Cable, 3 m (10 ft)	STANDARD STEM (FOR ELECTROSURGICAL PENCIL) 	SPLIT STEM (FOR FOOT CONTROL CORD USE) 
	0075	1

ELECTROSURGICAL PENCILS WITH E-Z CLEAN™



ELECTROSURGICAL PENCILS

Disposable Pencils with E-Z CLEAN™ Electrodes

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
		
Rocker Switch Pencil with E-Z CLEAN™ Electrosurgical Tip with 10' (3m) cord	0030	50
Rocker Switch Pencil with E-Z CLEAN™ Electrosurgical Tip with 15' (4.5m) cord	0030L	25
Rocker Switch Pencil with E-Z CLEAN™ Electrosurgical Tip and Holster with 10' (3m) cord	0030H	50
Rocker Switch Pencil with E-Z CLEAN™ Electrosurgical Tip and Holster with 15' (4.5m) cord	0030HL	25
		
Button Switch Pencil with E-Z CLEAN™ Electrosurgical Tip with 10' (3m) cord	0035	50
Button Switch Pencil with E-Z CLEAN™ Electrosurgical Tip with 15' (4.5m) cord	0035L	25
Button Switch Pencil with E-Z CLEAN™ Electrosurgical Tip and Holster with 10' (3m) cord	0035H	50
Button Switch Pencil with E-Z CLEAN™ Electrosurgical Tip and Holster with 15' (4.5m) cord	0035HL	25
		
Rocker Switch Pencil with Modified E-Z CLEAN™ Electrosurgical Tip, 15' (4.5m) cord	0036	50
Rocker Switch Pencil with Modified E-Z CLEAN™ Electrosurgical Tip and Holster, 15' (4.5m) cord	0036H	50
		
Button Switch Pencil with Modified E-Z CLEAN™ Electrosurgical Tip, 15' (4.5m) cord	0037	50
Button Switch Pencil with Modified E-Z CLEAN™ Electrosurgical Tip and Holster, 15' (4.5m) cord	0037H	50

MEGADYNE™ TELESCOPING SOFT TISSUE DISSECTORS



TELESCOPING SOFT TISSUE DISSECTORS

MEGADYNE™ Telescoping Soft Tissue Dissectors
With and without smoke evacuation

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 Telescoping Smoke Evacuation Soft Tissue Dissector	ME7251ST	6
 Telescoping Smoke Evacuation Soft Tissue Dissector, Modified	ME725M1ST	6
 Telescoping Soft Tissue Dissector	ME7251T	6
 Telescoping Soft Tissue Dissector, Modified	ME725M1T	6

REUSABLE PATIENT RETURN ELECTRODES



REUSABLE PATIENT RETURN ELECTRODES

MEGA SOFT™ Patient Return Electrodes

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 MEGA SOFT™ Universal Plus Patient Return Electrode 91 cm (36") L x 51 cm (20") W, 0.6 cm (1/4") thickness For patients 0.36 kg (0.8 lbs.) and up	0847	1
 MEGA SOFT™ Universal Plus Patient Return Electrode, Dual Cord 91 cm (36") L x 51 cm (20") W, 0.6 cm (1/4") thickness For patients 0.36 kg (0.8 lbs.) and up	0848	1
 MEGA SOFT™ Universal Patient Return Electrode 91 cm (36") L x 51 cm (20") W, 0.3 cm (1/8") thickness For patients 0.36 kg (0.8 lbs.) and up	0845	1
 MEGA SOFT™ Universal Patient Return Electrode, Dual Cord 91 cm (36") L x 51 cm (20") W, 0.3 cm (1/8") thickness For patients 0.36 kg (0.8 lbs.) and up	0846	1
 MEGA SOFT™ Patient Return Electrode 117 cm (46") L x 51 cm (20") W, 1.25 cm (1/2") thickness For patients 11.3 kg (25 lbs.) and up	0830	1
 MEGA SOFT™ Dual Cord Patient Return Electrode 117 cm (46") L x 51 cm (20") W, 1.25 cm (1/2") thickness For patients 11.3 kg (25 lbs.) and up	0835	1
 MEGA SOFT™ Pediatric Patient Return Electrode 66 cm (26") L x 30.5 cm (12") W, 1.3 cm (1/2") thickness For patients between 0.36 kg (0.8 lbs.) and 22.7 kg (50 lbs.)	0840	1
MEGA 2000™ Reusable Patient Return Electrode 91 cm (36") L x 51 cm (20") W For patients 11.3 kg (25 lbs.) and up	0800	1
MEGA 2000™ Sheath	0825	30

REUSABLE PATIENT RETURN ELECTRODES

MEGA SOFT™ Cables

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 MEGA SOFT™ Reusable DetachaCable™ Cord, Standard Connector, 2.4m (8.0 ft.)	M2K01	1
 MEGA SOFT™ Reusable DetachaCable™ Cord, Standard Connector, 4.4m (14.4 ft.)	M2K02	1
 MEGA SOFT™ Reusable DetachaCable™ Cord, Phone Plug Connector, 2.4m (8.0 ft.)	M2K03	1
 MEGA SOFT™ Reusable DetachaCable™ Cord, Phone Plug Connector, 4.4m (14.4 ft.)	M2K04	1
 MEGA SOFT™ Reusable DetachaCable™ Cord, Extended Phone Plug Connector, 4.4m (14.4 ft.)	M2K05	1
 MEGA SOFT™ Reusable DetachaCable™ Cord, Argon Beam Connector, 2.4m (8.0 ft.)	M2K06	1
Replacement Pigtail Cable	M2K07	1
 DetachaCable™ Cord, Compatibility, MEGA SOFT™ 2.4m (8.0 ft.)	M2K08	1
DetachaCable™ Cord, Compatibility, MEGA SOFT™ 4.4m (14.4 ft.)	M2K09	1
MEGA SOFT™ Patch Kit	PKIT001	1

DISPOSABLE PATIENT RETURN ELECTRODES



DISPOSABLE PATIENT RETURN ELECTRODES

Pre-Attached Cords

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 Adult, single plate, with 9' (2.7m) pre-attached cord, >33lbs (>15kg)	O850C	50
 Adult, dual plate, with 10' (3m) pre-attached cord, >33lbs (>15kg)	O855C	50
Adult, dual plate, with 15' (4.6m) pre-attached cord, >33lbs (>15kg)	O855CL	50
Adult, dual plate, with 9' (2.7m) pre-attached cord, >33 lbs (>15Kg)	O855CN	50
 Pediatric dual plate, with 3m (10') pre-attached cord, 6-33 lbs (2.7-15kg)	O865C	50

Reusable Cords

 Adult, dual plate, no cord, >33lbs (>15kg)	O855	50
---	------	----

SURGICAL SMOKE EVACUATORS



SURGICAL SMOKE EVACUATORS

MEGADYNE™ Smoke Evacuator and Accessories

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
		
MEGADYNE™ Smoke Evacuator	MESE1	1
		
MEGADYNE™ Filter	2550J	1
		
MEGADYNE™ Fluid Trap	2555J	10
		
MEGADYNE™ RF Sensor	2255J	1
		
MEGADYNE™ Connect Cable, 1m	2403J	1
		
MEGADYNE™ Connect Cable, 2.1m	2406J	1

SURGICAL SMOKE EVACUATORS

MINI VAC™ Smoke Evacuator and Accessories

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
		
MEGADYNE™ MINI VAC™ Smoke Evacuator, 220V	ECVV220	1
		
EZ Link RF Sensor	MGEZLINK01	1
		
MicroSafe Filter 1 pk	MGVS353	1
MicroSafe Filter 2 pk	MGVS35302	2
		
Fluid trap	MGVSFT10	10

SMOKE EVACUATION PENCILS



SMOKE EVACUATION PENCILS

MEGADYNE™ Telescoping Smoke Evacuation Pencil

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
		
MEGADYNE™ Telescoping Smoke Evacuation Pencil, 10-Foot (3m) Tubing	251010J	20
		
MEGADYNE™ Telescoping Smoke Evacuation Pencil, 15-Foot (4.5m) Tubing	251015J	20

SMOKE EVACUATION PENCILS

Smoke Evacuation Pencils and Disposables

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
		
ZIP PEN™ Ergonomic Smoke Evacuation Pencil, 10' (3m) Tubing	252510	20
ZIP PEN™ Ergonomic Smoke Evacuation Pencil, 10' (3m) Tubing with 22mm connector	252510EC	20
		
Nozzle 2.7"/6.8cm	2540J	10
		
Nozzle 5.2"/13.2cm	2560J	10
		
ULTRA VAC™ Telescopic Smoke Evacuation Pencil, 10' (3m) Tubing	211010J	20
ULTRA VAC™ Telescopic Smoke Evacuation Pencil, 10' (3m) Tubing with 22 mm connector	211010EC	20
		
Attacha Vac Smoke Evacuation Shroud	212009J	25
Attacha Vac Smoke Evacuation Shroud with 22mm Connector	212009EC	25

SMOKE EVACUATION PENCILS

Smoke Evacuation Accessories

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 MEGADYNE™ Laparoscopic Tubing	2390J	25
 MEGADYNE™ Speculum Tubing	2395J	25
 22mm Male to 10mm Male Connector	2145J	1
 Universal Connector	2150J	1

STAINLESS STEEL ELECTRODES & PENCILS



STAINLESS STEEL ELECTRODES & PENCILS

Stainless Steel Blade Electrodes

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 2.5 inch (6.35 cm), Blade Electrode	0312	12
 2.5 inch (6.35 cm), Modified Blade Electrode	0312M	12
 2.75 inch (7 cm), Blade Electrode	0312A	12
 2.75 inch (7 cm), Modified Blade Electrode	0312AM	12
 4 inch (10.2 cm), Blade Electrode	0314A	12
 6.5 inch (16.5 cm), Blade Electrode	0314	12
 6.5 inch (16.5 cm), Modified Blade Electrode	0314M	12

STAINLESS STEEL ELECTRODES & PENCILS

Stainless Steel Needle Electrodes

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 2.75 inch (7 cm), Needle Electrode	0313	12
 2.75 inch (7 cm), Modified Needle Electrode	0313M	12
 4 inch (10.2 cm), Modified Needle Electrode	0316AM	12
 6 inch (15.2 cm), Needle Electrode	0316	12
 6 inch (15.2 cm), Modified Needle Electrode	0316M	12

Stainless Steel Ball Electrodes

 5 inch (12.7 cm), Ball Electrode	0300	12
 2 inch (5.1 cm), Ball Electrode	0315	12

STAINLESS STEEL ELECTRODES & PENCILS

Disposable Pencils With Stainless Steel Electrodes

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
		
Rocker Switch Pencil with Stainless Steel Electrosurgical Tip	0038	50
Rocker Switch Pencil with Stainless Steel Electrosurgical Tip and Holster	0038H	50
		
Button Switch Pencil with Stainless Steel Electrosurgical Tip	0039	50
Button Switch Pencil with Stainless Steel Electrosurgical Tip and Holster	0039H	50
		
Foot Control Pencil with Stainless Steel Electrosurgical Tip and Holster	FC39H	50
		
Stainless Steel Tip Cleaner	TC100	50



PENCIL HOLSTER

SUCTION IRRIGATION DEVICES



SUCTION IRRIGATION DEVICES

Endopath™ Probe Plus™ II Handles

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
		
Probe Plus II Hand Controlled, Pistol Grip Handle	EPH02	6
		
Probe Plus II Hand Controlled, Pencil Grip Handle	EPH04	6

Endopath™ Probe Plus™ II Electrodes

		
Probe Plus II Hook Electrode, Hollow Tip (Shaft Diameter 5mm, Length 34 cm)	EPS01	6
Probe Plus II Hook Electrode, Hollow Tip (Shaft Diameter 5mm, Length 29 cm)	EPS05	6
		
Probe Plus II Spatula Electrode, Hollow Tip (Shaft Diameter 5mm, Length 34 cm)	EPS02	6
Probe Plus II Spatula Electrode, Hollow Tip (Shaft Diameter 5mm, Length 29 cm)	EPS06	6
		
Probe Plus II Right Angle Electrode, Hollow Tip (Shaft Diameter 5mm, Length 34 cm)	EPS03	6
Probe Plus II Right Angle Electrode, Hollow Tip (Shaft Diameter 5mm, Length 29 cm)	EPS07	6
		
Probe Plus II Curved Dissector Electrode, Hollow Tip (Shaft Diameter 5mm, Length 34 cm)	EPS04	6

SUCTION COAGULATORS




SUCTION COAGULATORS

Suction Coagulators

DESCRIPTION	FRENCH SIZE	CATALOG NUMBER	QTY PER BOX
			
Foot Control, 6" (15.2 cm) cannula, 10ft. (3.0 m) cord	10	004010J	10
Foot Control, 6" (15.2 cm) cannula, 10ft. (3.0 m) cord	10	004025	25

			
Hand Control, 6" (15.2 cm) cannula, 10ft. (3.0 m) cord	10	004125	25
Hand Control, 6" (15.2 cm) cannula, 10ft. (3.0 m) cord	8	004225	25
Hand Control, 6" (15.2 cm) cannula, 10ft. (3.0 m) cord	12	004325	25

LLETZ LOOP ELECTRODES




SINGLE USE LLETZ LOOP ELECTRODES

Square Loop

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
 10 mm Wide x 7 mm Deep with 11.8 cm Shaft Orange	0410	12

Round Loops

 10 mm Wide x 5 mm Deep with 11.8 cm Shaft Purple	0420	12
 15 mm Wide x 6 mm Deep with 11.8 cm Shaft Gray	0430	12
 20 mm Wide x 10 mm Deep with 11.8 cm Shaft Red	0440	12
 10 mm Wide x 10 mm Deep with 11.8 cm Shaft Yellow	0450	12
 15 mm Wide x 12 mm Deep with 11.8 cm Shaft Green	0460	12
 20 mm Wide x 12 mm Deep with 11.8 cm Shaft White	0470	12
 20 mm Wide x 15 mm Deep with 11.8 cm Shaft Blue	0480	12
 20 mm Wide x 8 mm Deep with 11.8 cm Shaft Tan	0490	12

POWERSTAR™

Bipolar Scissors



POWERSTAR™

POWERSTAR™

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
		
POWERSTAR™ 115mm TENOTOMY	BP100	1
POWERSTAR™ 175mm TENOTOMY	BP120	1
		
POWERSTAR™ FINE-TIP 140mm METZENBAUM	BP310	1
POWERSTAR™ FINE-TIP 175mm METZENBAUM	BP320	1
POWERSTAR™ FINE-TIP 230mm METZENBAUM	BP340	1
POWERSTAR™ FINE-TIP 280mm METZENBAUM	BP360	1
		
POWERSTAR™ 175mm METZENBAUM	BP520	1
POWERSTAR™ 230mm METZENBAUM	BP540	1
POWERSTAR™ 280mm METZENBAUM	BP560	1
		
POWERSTAR™ 170mm MAYO	BP720	1

POWERSTAR™

POWERSTAR™ Accessories

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
		
POWERSTAR™ CABLE 5m COAX 8mm	BP910	1
		
POWERSTAR™ CABLE 5m COAX 12mm	BP930	1
		
POWERSTAR™ CABLE 5m FIXED LEAD	BP940	1

ENDOPATH™

Bipolar Forceps



ENDOPATH™

ENDOPATH™

DESCRIPTION	CATALOG NUMBER	QTY PER BOX
		
Shaft Diameter 5mm, Jaw Macro	EBF01	1
Shaft Diameter 5mm, Jaw Micro	EBF02	1
Reusable coaxial pin connector bipolar cord, gray, 1 prong	EBC02	1
Active return cord for bipolar instruments fixed pin connector	EBC05	1

