



# Kepler Advanced

## Technical Specs

Loupe lens system	Kepler, but with higher light transmission than Kepler Kompakt
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System™	Yes (optional add-on)
Fine adjustable	Yes
Lockable Fine-focus available	No
Power up-system	No, not needed (come with 4 mags in one)
Magnification	Optional via click: 3.6x + 4.5x + 5.5x + 6.4x
Working distance	From 32 to 65 cm
Field of view (at 40 cm)	3.6x: 7 cm   4.5x: 6 cm   5.5x: 5 cm   6.4x: 4 cm
Depth of Field (at age 50 years)	3.6x: 8 cm   4.5x: 6 cm   5.5x: 5 cm   6.4x: 3.5 cm
Mount system	TTL
Weight in an Icon Medium frame	75,5 g
Loupe colour	Anthracite grey/Metallic grey
Loupe material	A carbon enhanced composite housing and a titanium ring for changing mag
Frames available	All