



Kepler Kompakt

Technical Specs

| | |
|----------------------------------|--|
| Loupe lens system | Kepler Kompakt |
| Lenses | Plano or Prescription |
| BlueLight Protection (BLP) | Yes (optional add-on) |
| Prisma Lens System™ | Yes (optional add-on) |
| Fine adjustable | Yes |
| Lockable Fine-focus available | Yes |
| Power up-system | Yes, front ocular is changeable |
| Magnification | 3.5x 4.6x 5.7x |
| Working distance | From 30 to 60 cm |
| Field of view (at 40 cm) | 3.5x: 8.5 cm 4.6x: 7.5 cm 5.7x: 6.5 cm |
| Depth of Field (at age 50 years) | 3.5x: 10 cm 4.6x: 8 cm 5.7x: 6 cm |
| Mount system | TTL |
| Weight in an Icon Medium frame | 75 g |
| Loupe colour | Grey Titanium housing. TTL black composite |
| Loupe material | Titanium and composite housing |
| Frames available | All (Except Carbon) |