



Kepler Kompakt

Technical Specifications

Variant	Kepler Kompakt v2 – 3.4x, 5.0x & 7.6x
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System	Yes (optional add-on)
Fine-focus available	Yes
Magnification	3.4x 5.0x 7.6x
Working distance	3.4x & 5.0x: From 33 cm to 66 cm 7.6x: From 39 cm to 66 cm
Field of View (at 45 cm)	3.4x: 9.5 cm 5.0x: 6.2 cm 7.6x: 4.3 cm
Depth of field (up to)	3.4x: 11.5 cm 5.0x: 8.0 cm 7.6x: 4.3 cm
Mount system	TTL (Through the lens)
Weight of two oculars with correction lenses	3.4x: 37.3 g 5.0x: 39.2 g 7.6x: 43.7 g
Loupe colour	Anthracite grey & Metallic grey
Loupe material	Carbon fibre reinforced composite and titanium
Frames available	All