



PowerGo

Technical Specifications

Control Unit	PowerGo
Material	Aluminium
Weight	155 g
Compatible LED Lights	Total Pure & Total Intense LED Lights
Connection type	USB-C
Light intensity settings	High, Medium, and low
Operating time	Varies with light intensity, 8-17 hours
Charging time	Approx. 2-2.5 hours
Battery Type	21,700 Lithium-ion rechargeable battery
Output rating	3VDC 0,5-2,0A
Input rating	5 VDC 2,0-2,5 A
Warranty	PowerGo Unit: 3 years PowerGo Battery: 1 year