

MatrixMIDFACE™ 1.3. Midface plating system

This description alone does not provide sufficient background for direct use of the product. Instruction in handling this product in accordance with the corresponding surgical technique (SE_807065) and by a surgeon experienced in handling this product is essential.


MatrixMIDFACE 1.3 Screws

Screw Type	Diameter (mm)	Lengths (mm)
Self-Tapping	1.3	3–8 10–12
Self-Drilling	1.3	3–6
Emergency	1.8	3–12




MatrixMIDFACE Drill Bits

Coupling	Diameter (mm)	Stop (mm)
J-Latch	1.0	4, 6, 8, 12
Mini Quick	1.0	4, 6, 8




Note: If a pilot hole is desired, use the appropriate 1.0 mm diameter drill bit for all screw lengths.

MatrixMIDFACE 1.3 Plates



Mini plates,
0.5 mm thick



MatrixMIDFACE 1.3 Modules & Sets



61.503.814 Module 1.3 MatrixMIDFACE, 3/3, with Lid, without Contents
01.503.816 1.3 MatrixMIDFACE Set, 3/3



61.503.601 Module 1.3 MatrixMIDFACE, 1/3, with Lid, without Contents, for use with sterile Implants
01.503.603 1.3 MatrixMIDFACE Set, 1/3 MatrixMIDFACE Set, 1/3, for use with sterile Implants


MatrixMIDFACE 1.3 vs. MatrixMIDFACE

Notes:

- MatrixMIDFACE 1.3 screws and MatrixMIDFACE dark blue mini plates are fully compatible with all MatrixMIDFACE instruments.
- MatrixMIDFACE 1.3 screws can only be used with dark blue MatrixMIDFACE mini plates.
- MatrixMIDFACE screws can be used with all MatrixMIDFACE plates, except for dark blue MatrixMIDFACE mini plates.
- MatrixMIDFACE dark blue mini plates have a smaller hole-to-hole distance compared to comparable MatrixMIDFACE light blue mini plates.

System	Screw	Plate				
		Micro Silver 0.4 mm	Mini Light Blue 0.5 mm	Mini Dark Blue 0.5 mm	Medium Pink 0.7 mm	Large Gold 0.8 mm
MatrixMIDFACE 1.3	Self-Tapping 1.3					
	Self-Drilling 1.3			X		
	Emergency 1.7					
MatrixMIDFACE	Self-Tapping 1.5					
	Self-Drilling 1.5	X	X		X	X
	Emergency 1.8					

MatrixMIDFACE 1.3



MatrixMIDFACE

