

BECKER
► **ICeline**



MOBILE FURNITURE

EASY TO MOVE, PRACTICAL AND STURDY

05/2024

BECKER ICEline

MOBILE FURNITURE BY BECKER!

Good hospital furniture is characterised by its suitability for everyday use. The new BECKER ICEline offers a wide range of tried and tested mobile assistants for all sectors. These are all fully and suitably equipped for their intended use and make work a lot easier.

Using BECKER hook trolleys, your clinical equipment is always in the right place. Classical Becker Equipment will accompany your everyday clinical work reliably for many years. Our mobile BECKER sterile goods components are ideal for the transportation, supply and disposal of sterile and soiled items.

Furthermore, BECKER's selection of instrument and service trolleys provides an important link in the logistics chain of the entire medical and hospital sector. These components can be combined freely, are designed as a modular system and can be easily configured for specific purposes from a number of basic models.

To put it simply: The main objective of the BECKER ICEline is to streamline and optimise OP, outpatient and CSSD processes while ensuring hygiene conditions at all times. For you as a user, this means greater flexibility and convenience in all fields of application.

ICE. ICE. Baby.



CONTENTS

- ▶ 04. Unbeatable advantages
- ▶ 06. Hook trolleys
- ▶ 14. Equipment
- ▶ 26. Container trolleys
- ▶ 36. Instrument trolleys
- ▶ 42. ISO-trolleys
- ▶ 48. Service trolleys



UNBEATABLE ADVANTAGES

WHY BECKER ICELINE?

► IT IS STRONG, ROBUST AND MODULAR

Ideal prerequisites for effective work assistance.

► NEW DESIGN

Brings ease, functionality and modern design to your hospital.

► TOP QUALITY AND CLEVER CONSTRUCTION

Here we stick to our principles: optimally designed and with a sturdy construction to ensure a long service life.

► WASHABLE



All models with this marking are suitable for use in machine washers.

► STEEL

We use the material we trust to ensure hygiene and a long service life. Our signal for the future.



CNS: AN ALLROUND TALENT

OUR FOCUS IS ON SUSTAINABILITY!

CHROM NICKEL STEEL

- ▶ pleasant appearance and ideal workability
- ▶ durable, ageless and reusable



SUSTAINABILITY IS IMPORTANT TO US!


By using steel, we aim to create durable, robust and economical products without having to constantly create new models.

What's more, thanks to its unique recycling properties, steel can be returned to the material cycle almost indefinitely.

WHAT EXACTLY IS CNS 1.4301?

This is the most commonly-used form of austenitic chrom nickel steel alloys. Its outstanding properties are:

- ▶ excellent durability
- ▶ corrosion resistance thanks to chromium component
- ▶ withstands temperatures up to 600 °C



HOOK TROLLEYS

QUALITY AND TECHNICAL DETAILS

OUR ARTICLE NUMBER: TH _ _ _ _ _

Tube

Hooks

Basic model TH00



Can be fitted with sterile goods baskets, shoe holders and suture material shelves

Basic model TH03



Can be fitted with wire baskets

QUALITY

Made of square chrom nickel steel tubing, 25 x 25 mm (models TH00 and TH03), 40 x 40 mm (model TH05), material no. 1.4301 (AISI 304) with a 320 grit surface finish.

HOOKS

To hang accessories onto the frame. Can be flexibly combined and will be in stock for many years.

Basic model TH05



Can be fitted with
M9 and M11 baskets

WASHABLE

Hook-rack trolleys with these markings are suitable for machine washers.



CASTORS

All castors are of premium design: electrically conductive, non-marking and with thread guard.

lockable polymer
double castors Ø 65 mm



lockable, washable polymer
swivel castors Ø 125 mm



HOOK TROLLEYS TH00

For sterile item baskets and other accessories



Hook trolleys, single-sided

3 pairs of hooks


Article number
TH00030

4 lockable polymer
double castors Ø 65 mm,
wall protectors

 645 x 540 x 1450 mm

Article number
TH00037

4 lockable, washable polymer
swivel castors Ø 125 mm,
wall protectors


 645 x 540 x 1535 mm



4 pairs of hooks

Article number
TH00040

4 lockable polymer
double castors Ø 65 mm,
wall protectors

 645 x 540 x 1830 mm

Article number
TH00047

4 lockable, washable polymer
swivel castors Ø 125 mm,
wall protectors

 645 x 540 x 1915 mm




Hook trolleys, double-sided

3 pairs of hooks


Article number
TH00060

4 lockable polymer
double castors Ø 65 mm,
wall protectors

 645 x 710 x 1450 mm

Article number
TH00067

4 lockable, washable polymer
swivel castors Ø 65 mm,
wall protectors


 645 x 710 x 1535 mm



4 pairs of hooks


Article number
TH00080

4 lockable polymer
double castors Ø 65 mm,
wall protectors

 645 x 710 x 1830 mm

Article number
TH00087

4 lockable, washable polymer
swivel castors Ø 125 mm,
wall protectors

 645 x 710 x 1915 mm



All dimensions  stated as width x depth x height (W/D/H)

	Description	Article number	 W/D/H
	Sterile goods basket 1/2-StU	120100	575 x 280 x 135 mm
	Wire partition for sterile goods basket 1/2-StU	120103	280 x 120 mm
	Sterile goods basket 1/2-StU with slanted sides and lower front access	120110	575 x 280 x 135 mm
	Wire partition for sterile goods basket 1/2-StU with slanted sides and lower front access	120113	255 x 120 mm
	Label plate for Sterile goods basket 1/2-StU	120116	100 x 80 mm
	Sterile goods basket 1-StU	120120	575 x 280 x 265 mm
	Wire partition for sterile goods basket 1-StU with or without lower front access	120123	280 x 240 mm
	Sterile goods basket 1-StU with lower front access 300 x 120 mm	120130	575 x 280 x 265 mm
	Label plate for Sterile goods basket 1-StU	120146	100 x 100 mm
	Pair of slanted holders for mounting a sterile goods basket or a shoe rack at an angle	120160	35 x 50 x 115 mm
	Shoe rack for storing 6 pairs of OR shoes	120430	600 x 265 x 190 mm
	Suture material shelf for 20 normal suture material boxes	120530	600 x 165 x 290 mm



Abb. Slanted sterile goods basket

Note on sterile goods baskets:

Centred wire loops, welded top and bottom frames, bottom frame designed for stacking, tested by TÜV (German technical inspection association), design conforms with DIN 58952.

Load capacity: 1/2-StU basket approx 20 kg. 1-StU basket approx 30 kg.

HOOK TROLLEYS TH00

Examples of configurations



Description

Hook trolley
double-sided with 4 pairs of hooks incl.

- ▶ steril goods basket ½-StU with slanted access and wire partition
- ▶ ½-StU with wire partition and labeling plate
- ▶ 1-StU with wire partition
- ▶ 1-StU with opening and labeling plate

4 lockable, washable polymer swivel castors Ø 125 mm, wall protectors

Article number

TH00987

W/D/H

645 x 710 x 1915 mm



Description

OR shoe rack with
3 shoe holders incl.

- ▶ 3 inclined adjuster
- ▶ 3 shoe holders for holding 6 pairs of surgical shoes each (18 pairs in total)

4 lockable, washable polymer swivel castors Ø 125 mm, wall protectors

Article number

TH00137

W/D/H

645 x 540 x 1535 mm



Description

OR shoe rack with
4 shoe holders incl.

- ▶ 4 inclined adjuster
- ▶ 4 shoe holders for holding 6 pairs of surgical shoes each (24 pairs in total)

4 lockable, washable polymer swivel castors Ø 125 mm, wall protectors

Article number

TH00147

W/D/H

645 x 540 x 1915 mm



Description

Hook trolley with
4 pairs of hooks incl.

- ▶ 2 suture material cabinets

4 lockable, washable polymer double castors Ø 65 mm, wall protectors

Article number

TH00944

W/D/H

645 x 540 x 1830 mm

A tall, narrow metal wire shelving unit with five shelves and four casters. The unit has a simple, functional design with a dark metal frame and wire mesh shelves. It is positioned on a light-colored surface against a plain background.

Accessories TH03

A three-tier wire shelving unit on wheels, model TH05127, is shown. The unit has a silver-colored metal frame and wire mesh shelves. It is equipped with four casters, two of which are swivel casters with red wheels. A red arrow points to the unit, and the text "Can be fitted" is displayed above it.

TH05127

11

”

*The magic
is in the
detail.*

Theodor Fontane, German writer


HOOK TROLLEYS TH05

Can be fitted with M9 and M11 baskets






TH05137

illustration with
2 additional wire divider
Article number.: 410

Description	Article number	 W/D/H
OR trolley with 3 pairs of hooks incl.		
<ul style="list-style-type: none"> ▶ 3 M9 baskets 	TH05137	845 x 660 x 1730 mm
<ul style="list-style-type: none"> or ▶ 3 M11 baskets 	TH05237	1045 x 660 x 1730 mm
4 lockable, washable polymer swivel castors Ø 125 mm		

Accessories TH05

	Description	Article number	 W/D/H
	M9 basket	410 67 090	845 x 560 x 180/400 mm
	M11 basket	410 67 110	1045 x 560 x 180/400 mm
	wire partition	410	520/335 x 360/150 mm

EQUIPMENT

QUALITY AND TECHNICAL DETAILS

OUR ARTICLE NUMBER: TE _ _ _ _ _

Tube

Equipment



Mayo tables
TE01



Arm rest TE02



Kick buckets TE03



Bowl stand TE04 and
instrument drop TE05

14.

QUALITY

Made of chrom nickel steel tubes and sheets, material no. 1.4301 (AISI 304) with a 320 grit surface finish.

WASHABLE

Equipment with these markings are suitable for machine washers.




Surgeon's step
TE06

CASTORS

All castors are of premium design: electrically conductive, non-marking and with thread guard.

lockable polymer double castors, Ø 65 mm



lockable polymer double castors, Ø 100 mm



lockable, washable polymer swivel castors, Ø 125 mm



FEET

Antistatic conductive feet



Infusion stands
TE07



Linen and waste
trolleys TE08



Sterilization wrap
stands TE09

MAYO TABLES TE01



TE01100

Description	Article number	W/D/H
manually height adjustable 835 - 1265 mm		
Deep-drawn plate with rounded corners		
▶ plate size 600 x 400 mm	TE01100	680 x 500 x 840 mm
or		
▶ plate size 700 x 500 mm	TE01110	700 x 500 x 840 mm
4 lockable, polymer double castors Ø 65 mm, wall protectors		



TE01220

Description	Article number	W/D/H
oil-hydraulically height adjustable 920 - 1305 mm		
Deep-drawn plate with rounded corners, 360° rotatable and lockable		
▶ plate size 600 x 400 mm	TE01200	700 x 490 x 925 mm
or		
▶ plate size 700 x 500 mm	TE01210	700 x 500 x 925 mm
or		
▶ plate size 850 x 600 mm	TE01220	850 x 600 x 925 mm
3 lockable, polymer double castors Ø 65 mm		

Height adjustment range is changed due to larger castors




TE01214

oil-hydraulically height adjustable 960 - 1345 mm		
Deep-drawn plate with rounded corners, 360° rotatable and lockable		
▶ plate size 600 x 400 mm	TE01204	700 x 490 x 965 mm
or		
▶ plate size 700 x 500 mm	TE01214	700 x 500 x 965 mm
or		
▶ plate size 850 x 600 mm	TE01224	850 x 600 x 965 mm
3 lockable, polymer double castors Ø 100 mm		

MAYO TABLES TE01




Description	Article number	 W/D/H
<p>oil-hydraulically height adjustable 920 - 1305 mm</p> <p>Deep-drawn plate with rounded corners, can be swiveled 360° and locked in all positions</p> <p>▶ plate size 1220 x 620 mm</p> <p>4 lockable, polymer double castors Ø 65 mm</p>	TE01300	1220 x 620 x 925 mm

Height adjustment range is changed due to larger castors



<p>oil-hydraulically height adjustable 960 - 1345 mm</p> <p>Deep-drawn plate with rounded corners, can be swiveled 360° and locked in all positions</p> <p>▶ plate size 1220 x 620 mm</p> <p>4 lockable, polymer double castors Ø 100 mm</p>	TE01304	1220 x 620 x 965 mm
--	---------	------------------------



Description	Artikel-Nr	 W/D/H
<p>oil-hydraulically height adjustable 920 - 1305 mm</p> <p>Deep-drawn plate with rounded corners, top tray can be swiveled 360° and locked in all positions</p> <p>▶ plate size top and bottom 1220 x 620 mm</p> <p>4 lockable, polymer double castors Ø 65 mm</p>	TE01310	1220 x 620 x 925 mm

Height adjustment range is changed due to larger castors



<p>oil-hydraulically height adjustable 960 - 1345 mm</p> <p>Deep-drawn plate with rounded corners, top tray can be swiveled 360° and locked in all positions</p> <p>▶ plate size top and bottom 1220 x 620 mm</p> <p>4 lockable, polymer double castors Ø 100 mm</p>	TE01314	1220 x 620 x 965 mm
--	---------	------------------------

ARMREST TE02



Description

Arm rest
manually height adjustable
790 - 1220 mm

- ▶ without cushion
- ▶ with adjustable table top
450 x 325 mm

5 lockable, polymer double castors
Ø 65 mm

Article number

TE02100

Ø/H

Ø600,
790 mm



components

(please order separately)

Cushion for arm rest (black)

- ▶ electrically conductive

TE0210K

W/D/H

450 x 325 x
40 mm



Description

Hand and arm rest
manually height adjustable
675 - 975 mm

- ▶ without cushion
- ▶ with radiolucent table top
470 x 450 mm, adjustable in all
directions thanks to ball-type
swivel joint

4 lockable, polymer double castors
Ø 65 mm, wall protectors

Article number

TE02200

W/D/H

500 x 540 x
675 mm



components

(please order separately)

Cushion for hand and arm rest (black)


- ▶ radiolucent
- ▶ electrically conductive

TE0220K


470 x 450 x
40 mm

KICK BUCKETS TE03




Description	Article number	 Ø/H
<ul style="list-style-type: none"> with bucket 15 l, chrome nickel steel, sterilizable <p>4 lockable, polymer double castors Ø 65 mm</p>	TE03100	Ø435, 450 mm




Description	Article number	 Ø/H
<ul style="list-style-type: none"> with bucket 10 l, plastic, white <p>4 lockable, polymer double castors Ø 65 mm</p>	TE03120	Ø435, 370 mm




Description	Article number	 W/D/H
<ul style="list-style-type: none"> with bucket 15 l, chrome nickel steel, sterilizable with push handle <p>4 lockable, polymer double castors Ø 65 mm</p>	TE03200	435 x 650 x 965 mm



Description	Article number	 W/D/H
<ul style="list-style-type: none"> with bucket 10 l, plastic, white with push handle <p>4 lockable, polymer double castors Ø 65 mm</p>	TE03220	435 x 650 x 965 mm




Description	Article number	 W/D/H
<p>manually height adjustable 430 - 920 mm</p> <ul style="list-style-type: none"> without bucket for holding a 15 liter chromium-nickel steel bucket with push handle <p>4 lockable, polymer double castors Ø 65 mm</p>	TE03300	510 x 510 x 1010 mm

Accessories TE03


	Description	Article number
	plastic bucket 10 l, white	130129
	bucket 15 l, chrome nickel steel, sterilizable	130128
	lid for 15 l chrome nickel steel bucket, sterilizable	130127

BOWL-HOLDER TE04





Description	Article number	 W/D/H
<ul style="list-style-type: none"> ▶ with chrome nickel steel bowl 8 l ▶ push handle <p>4 lockable, polymer double castors Ø 65 mm</p>	TE04100	550 x 560 x 980 mm




Description	Article number	 W/D/H
<ul style="list-style-type: none"> ▶ with 2 chrome nickel steel bowls 8 l ▶ push handle <p>4 lockable, polymer double castors Ø 65 mm</p>	TE04120	500 x 790 x 980 mm

Accessories TE04


	Description	Article number	 Ø/H
	chrome nickel steel bowl 8 l	130148	Ø220/360, 150 mm

INSTRUMENT DROP TE05





Description	Article number	 W/D/H
<p>Linen hamper with linen bag</p> <ul style="list-style-type: none"> ▶ ring Ø 400 mm, chrome nickel steel ▶ push handle <p>4 lockable, polymer double castors Ø 65 mm</p>	TE05100	550 x 560 x 980 mm



Description	Article number	 W/D/H
<p>Instrument receptacle with shelf 518 x 318 mm</p> <ul style="list-style-type: none"> ▶ without wire tray ▶ push handle <p>4 lockable, polymer double castors Ø 65 mm</p>	TE05200	550 x 560 x 980 mm

Accessories TE05

	Description	Article number	 W/D/H
	<p>wire tray, chrome nickel steel with drop handle</p>	130178	480 x 250 x 50 mm

SURGEON`S STEPS TE06

surgeon`s single step



Description

precisely fitting, removable and non-slip plate

▶ Loadable up to 150 kg

Article number

TE0610

W/D/H

600 x 350 x 130 mm



TE0612

600 x 350 x 220 mm



TE0614

600 x 350 x 300 mm

surgeon`s double step



Description

precisely fitting, removable and non-slip plate

▶ Loadable up to 150 kg

Article number


TE0622

W/D/H

600 x 650 x 220/440 mm

SURGEON'S STEPS TE06




Description	Article number	 W/D/H
Trolley to store 4 surgeon's single steps (without surgeon's steps)		
<p>mobil</p> <p>4 lockable, washable polymer swivel castors Ø 125 mm</p> <p>or</p> <p>wall-mounted system</p>	TE06647	620 x 420 x 1600 mm
	TE0654	600 x 430 x 1450 mm

TE06647




INFUSION STANDS TE07




Description	Article number	 Ø/H
<p>manually height adjustable</p> <p>1360 - 2440 mm</p> <ul style="list-style-type: none"> ▶ with 4 hooks ▶ PP drip catcher <p>5 lockable, polymer double castors Ø 65 mm</p>	TE07100	Ø660, 1360 mm




Description	Article number	 Ø/H
<p>manually height adjustable</p> <p>1360 - 2440 mm</p> <ul style="list-style-type: none"> ▶ with 4 hooks ▶ push handle <p>5 lockable, polymer double castors Ø 65 mm</p>	TE07200	Ø660, 1360 mm


LINEN AND WASTE TROLLEYS TE08



Description	Article number	 W/D/H
<ul style="list-style-type: none"> ▶ without bags ▶ for holding 1 bag ▶ with clipring and push handle ▶ removable bottom plate <p>4 lockable, polymer double castors Ø 65 mm</p>	TE08100	430 x 430 x 900 mm



Description	Article number	 W/D/H
<ul style="list-style-type: none"> ▶ without bags ▶ for holding 2 bag ▶ with 2 cliprings and push handle ▶ removable bottom plate <p>4 lockable, polymer double castors Ø 65 mm</p>	TE08200	970 x 430 x 900 mm
<div>components</div> <p>(please order separately)</p> <ul style="list-style-type: none"> ▶ Lid for linen and waste trolley, chrome nickel steel 	TE0880K	370 x 370 x 70 mm
<div>spare part</div> <ul style="list-style-type: none"> ▶ Clipring with push handle, chrome nickel steel 	TE0890K	355 x 390 x 35 mm



STERILIZATION WRAP STANDS TE09



Description

- ▶ with 4 storage frames
- ▶ for wrap sizes up to max. 1500 x 2050 mm

4 lockable, polymer double castors
Ø 65 mm, wall protectors

Article number

TE09400

W/D/H

1600 x 595 x
775/1220 mm



components

(please order separately)

- ▶ Bottom plate
for hanging in sterilization wrap
stands, with 4 storage frames,
chrome nickel steel

TE0940K

1370 x 460 x
20 mm



Description

- ▶ with 5 storage frames
- ▶ for wrap sizes up to max. 1500 x 2050 mm

4 lockable, polymer double castors
Ø 65 mm, wall protectors

Article number

TE09500

W/D/H

1600 x 760 x
775/1220 mm



components

(please order separately)

- ▶ Bottom plate
for hanging in sterilization wrap
stands, with 5 storage frames,
chrome nickel steel

TE0950K

1370 x 620 x
20 mm

CONTAINER TROLLEYS

QUALITY AND TECHNICAL DETAILS

OUR ARTICLE NUMBER: TC _ _ _ _ _

Tube

Container



Sterile goods
transportation trolleys TC07



Sterile goods
shelf trolleys TC03



Sterile goods
container trolleys TC01

26.

QUALITY

Made of chrom nickel steel tubes and sheets, material no. 1.4301 (AISI 304) with a 320 grit surface finish.

BASIC MODEL TC07

Sturdy base frame with cabinet: double-walled side panels, top and a pair of lockable swinging doors opening to an angle of 270 degrees and which can be locked in that position; two push handles and wall-protection bumpers; to be fitted with $\frac{3}{4}$ -extension wire shelves or height-adjustable single-shelled CNS shelves.

WASHABLE

Container trolleys with these markings are suitable for machine washers.



CASTORS

All castors are of premium design: electrically, conductive, non-marking and with thread guard.

lockable, polymer double castors, Ø 100 mm



lockable, washable polymer swivel castors, Ø 125 mm



lockable, swivel castors, Ø 125 mm



lockable, washable swivel castors, Ø 160 mm, chrome nickel steel



**Sterile goods
shelf/case cart trolleys TC05**



**Sterile goods
preparation trolleys TC09**

STERILE GOODS CONTAINER TROLLEYS TC01



Description

for sterile goods container
600/300 mm

- ▶ with 6 wire trays

Article number

TC01107

W/D/H

645 x 420 x
1530 mm

or





- ▶ with each 6 wire trays
on side
(12 trays in total)

4 lockable, washable polymer swivel
castors Ø 125 mm

TC01207

645 x 770 x
1530 mm


Accessories TC01

	Description	Article number	 W/D/H
	wire tray, chrome nickel steel	130478	620 x 320 x 175 mm


 W/D/H

STERILE GOODS SHELF TROLLEYS TC03



Description	Article number	 W/D/H
Clear shelf space 390 mm		
▶ with 4 fix welded shelves	TC03407	1340 x 730 x 1520 mm
or		
▶ with 4 fix welded and perforated shelves	TC03417	1340 x 730 x 1520 mm
4 lockable, washable polymer swivel castors Ø 125 mm, wall protectors		



Description	Article number	 W/D/H
Clear shelf space 285 mm		
▶ with 5 fix welded shelves	TC03507	1340 x 730 x 1520 mm
or		
▶ with 5 fix welded and perforated shelves	TC03517	1340 x 730 x 1520 mm
4 lockable, washable polymer swivel castors Ø 125 mm, wall protectors		

STERILE GOODS SHELF TROLLEYS / CASE CARTS TC05



Description

with 10 wire storage trays
chrome nickel steel, for holding

- ▶ 8 x ½-StU
- ▶ 2 x 1-StU
- ▶ max. 12 x ½-StU
(2 additional chrome nickel steel wire trays required)

4 lockable, washable polymer swivel castors Ø 125 mm, wall protectors

Article number

TC05107

W/D/H

820 x 690 x 1360 mm



Description

with 15 wire storage trays
chrome nickel steel, for holding

- ▶ 12 x ½-StU
- ▶ 3 x 1-StU
- ▶ max. 18 x ½-StU
(3 additional chrome nickel steel wire trays required)

4 lockable, washable polymer swivel castors Ø 125 mm, wall protectors



Article number

TC05157

W/D/H

1180 x 690 x 1360 mm

Accessories TC05

	Description	Article number	W/D/H
	Wire storage tray in StU measure, chrome nickel steel	160020	326 x 600 x 56 mm
	Identification plate can be hooked into wire storage tray, chrome nickel steel	160030	58 mm lang

STERILE GOODS SHELF TROLLEYS / CASE CARTS TC05



Description

with 4 shelves
top: with raised edge,
on tube frame
center: 2 shelves with raised edge
for hanging in bolts
bottom: with raised edge,
welded in place

- ▶ for holding OP accessories
- ▶ and/or 4 x ½-StU and 2 x 1-StU
- ▶ max. 6 x ½-StU and 2 x 1-StU

4 lockable, washable polymer swivel castors Ø 125 mm, wall protectors

Article number

W/D/H

TC05547

825 x 725 x
1280 mm

components

- ▶ Perforated flat shelf for case cart, for hanging in bolt, chrome nickel steel

TC0561K

825 x 725 x
25 mm

- ▶ Shelf on front and rear side with raised edge for case cart, usable on both sides, for hanging in bolt, chrome nickel steel

TC0563K

825 x 725 x
40 mm

- ▶ Shelf with raised edge for case cart, for hanging in bolt, chrome nickel steel

TC0564K

825 x 725 x
30 mm

- ▶ Push handle for hanging in bolt, chrome nickel steel

TC0579K

583 x 105 x
173 mm
(Trolley becomes
80 mm wider)

- ▶ MAGNES side panel for holding magnetic accessories, for hanging in bolt, magnetic chrome nickel steel

TC0570K

583 x 25 x
173 mm

- ▶ MAGNES side panel with standard rail for hanging in bolt, magnetic chrome nickel steel

TC0574K

583 x 75 x
173 mm
(Trolley becomes
50 mm wider)

- ▶ Panel with hooks for holding sterilisation baskets for hanging in bolt, chrome nickel steel

TC0578K

583 x 45 x
173 mm
(Trolley becomes
20 mm wider)

STERILE GOODS TRANSPORTATION TROLLEYS TC07



Description

with 10 wire storage trays (¼ pull-out)
for holding

- ▶ 8 x ½-StU and
- ▶ 2 x 1-StU
- ▶ max. 12 x ½-StU
(2 additional chrome nickel steel
wire shelves ¼ pull-out required)
- ▶ 2 lockable swinging doors
- ▶ 2 push handles

4 lockable swivel castors
Ø 125 mm

or

4 lockable, washable swivel castors
Ø 160 mm, chrome nickel steel



Article number

W/D/H

TC07306

850 x 720 x
1460 mm

TC07309

850 x 720 x
1500 mm



Description

with 15 wire storage trays (¼ pull-out)
for holding

- ▶ 12 x ½-StU and
- ▶ 3 x 1-StU
- ▶ max. 18 x ½-StU
(3 additional chrome nickel steel
wire shelves ¼ Pull-out required)
- ▶ 2 lockable swinging doors
- ▶ 2 push handles

4 lockable swivel castors
Ø 125 mm

or

4 lockable, washable swivel castors
Ø 160 mm, chrome nickel steel



Article number

W/D/H




TC07356

1190 x 720 x
1460 mm

TC07359

1190 x 720 x
1500 mm



	Description	Article number	 W/D/H
	Wire shelf 1/4 pull-out, chrome nickel steel	160010	326 x 600 x 56 mm
	Identification plate can be hooked into wire storage tray, chrome nickel steel	160030	58 mm length

Special feature of the 1/4 pull-out wire shelf

A special proprietary design feature makes this safe and practical to use

The sides of the 1/4 pull-out wire shelf form rails which slide over bolts in the trolley frame.

Asymmetrically located counteracting rails reliably prevent the extended shelf from tipping when a (heavy) load is placed on it.

To store a sterile goods container in the trolley, the wire shelf is pulled out as far as the 1/4 pull-out will allow, the end of the container is rested on the end of the shelf and is then pushed into the rack, also moving the wire shelf moves back into its original position.

Conversely, in order to remove a container, it must be pulled out towards the user. The wire shelf moves outwards with the container and provides a support so that the container can be lifted by the handles.

Simple and ergonomical !



STERILE GOODS TRANSPORTATION TROLLEYS TC07



Description

with 2 single-walled and removable shelves for holding

- ▶ max. 6 x ½-StU
- ▶ 2 lockable swinging doors
- ▶ 2 push handles

4 lockable swivel castors
Ø 125 mm

or

4 lockable, washable swivel castors
Ø 160 mm, chrome nickel steel



Article number

W/D/H

TC07506

850 x 720 x
1460 mm

TC07509

850 x 720 x
1500 mm



components

- ▶ Shelf for transport trolley for sterile goods, single walled and removable, chrome nickel steel

TC0750K

850 x 720 x
20 mm



Description

with 2 single-walled and removable shelves for holding

- ▶ max. 9 x ½-StU
- ▶ 2 lockable swinging doors
- ▶ 2 push handles

4 lockable swivel castors
Ø 125 mm

or

4 lockable, washable swivel castors
Ø 160 mm, chrome nickel steel



Article number

W/D/H

TC07556

1190 x 720 x
1460 mm

TC07559

1190 x 720 x
1500 mm



components


- ▶ Shelf for transport trolley for sterile goods, single walled and removable, chrome nickel steel

TC0755K


1190 x 720 x
20 mm

STERILE GOODS PREPARATION TROLLEYS TC09




Description	Article number	 W/D/H
with table top <ul style="list-style-type: none"> ▶ table top format 600 x 380 x 30 mm ▶ wire storage tray, chrome nickel steel flat welded shelf at the bottom 4 lockable, polymer double castors Ø 100 mm, wall protectors	TC09104	685 x 465 x 900 mm






Description	Article number	 W/D/H
with table top <ul style="list-style-type: none"> ▶ table top format 605 x 640 x 30 mm ▶ 2 wire storage tray flat welded shelf at the bottom 4 lockable, polymer double castors Ø 100 mm, wall protectors	TC09124	685 x 645 x 900 mm



Description	Article number	 W/D/H
with table top and raised edge <ul style="list-style-type: none"> ▶ table top format 665 x 380 x 30 mm ▶ 2 wire storage tray flat welded shelf at the bottom 4 lockable, polymer double castors Ø 100 mm, wall protectors	TC09144	685 x 465 x 900 mm

Accessories TC09

	Description	Article number	 W/D/H
	Wire storage tray in StU measure, chrome nickel steel	160020	326 x 600 x 56 mm
	Identification plate can be hooked into wire storage tray,	160030	58 mm length

INSTRUMENT TROLLEYS

QUALITY AND TECHNICAL DETAILS

OUR ARTICLE NUMBER: TI _ _ _ _ _

Tube

Instrument

36.



3 types of shelves:



Flat shelf

shelf with raised edge

shelf with gallery



welded in place



removable

QUALITY

Made of chrom nickel steel tubes and sheets, material no. 1.4301 (AISI 304) with a 320 grit surface finish.

WASHABLE

Instrument trolleys with these markings are suitable for machine washers.



CASTORS

All castors are of premium design: electrically conductive, non-marking and with thread guard.


lockable polymer double castors, Ø 65 mm




lockable polymer swivel castors, Ø 75 mm



lockable polymer double castors, Ø 100 mm



lockable, washable polymer swivel castors, Ø 125 mm



H. w. castors 75, 100	650	800	900	1000
H. w. castors 65	610	760	860	960
H. w. castors 125	685	835	935	1035
h₂	390	540	640	740
h₃	180	255	305	355
h₄ (with drawer)	245	395	495	595

B, T	450	650	750	950	1250
b, t	390	590	690	890	1190
GB, GT	530	730	830	1030	1330
GBS		805	905	1105	1405

INSTRUMENT TROLLEYS

All the models shown on this page are equipped with double castors Ø 100 mm

▶ top shelf: flat, welded on
bottom shelf: flat, welded on



size ▶	650 x 450 mm	750 x 650 mm	950 x 650 mm	1250 x 650 mm
height ▼	Article number			
650 mm	TI112 10104	TI312 10104	TI512 10104	TI712 10104
800 mm	TI132 10104	TI332 10104	TI532 10104	TI732 10104
900 mm	TI152 10104	TI352 10104	TI552 10104	TI752 10104
1000 mm	TI172 10104	TI372 10104	TI572 10104	TI772 10104

▶ top shelf: with raised edge, on tubular frame
bottom shelf: flat, welded on



size ▶	650 x 450 mm	750 x 650 mm	950 x 650 mm	1250 x 650 mm
height ▼	Article number			
650 mm	TI112 R0104	TI312 R0104	TI512 R0104	TI712 R0104
800 mm	TI132 R0104	TI332 R0104	TI532 R0104	TI732 R0104
900 mm	TI152 R0104	TI352 R0104	TI552 R0104	TI752 R0104
1000 mm	TI172 R0104	TI372 R0104	TI572 R0104	TI772 R0104

▶ top shelf: with gallery, on tubular frame
bottom shelf: flat, welded on



size ▶	650 x 450 mm	750 x 650 mm	950 x 650 mm	1250 x 650 mm
height ▼	Article number			
650 mm	TI112 G0104	TI312 G0104	TI512 G0104	TI712 G0104
800 mm	TI132 G0104	TI332 G0104	TI532 G0104	TI732 G0104
900 mm	TI152 G0104	TI352 G0104	TI552 G0104	TI752 G0104
1000 mm	TI172 G0104	TI372 G0104	TI572 G0104	TI772 G0104

Load-bearing capacity per shelf
for table size 650x450, 750x650, 950x650

- ▶ 60 kg per shelf plate (when 2 plates + rollers Ø 100 mm)
- ▶ 90 kg per shelf plate (when 2 plates + rollers Ø 125 mm)

Load-bearing capacity per shelf
for table size 1250x650

- ▶ 50 kg per shelf plate (when 2 plates + rollers Ø 100 mm)
- ▶ 50 kg per shelf plate (when 2 plates + rollers Ø 125 mm)

INSTRUMENT TROLLEYS

with double castors Ø 100 mm


- ▶ top shelf: with raised edge, on tubular frame
- bottom shelf: with raised edge, on tubular frame



size ▶	650 x 450 mm	750 x 650 mm	950 x 650 mm	1250 x 650 mm
height ▼	Article number			
650 mm	TI112 R0R04	TI312 R0R04	TI512 R0R04	TI712 R0R04
800 mm	TI132 R0R04	TI332 R0R04	TI532 R0R04	TI732 R0R04
900 mm	TI152 R0R04	TI352 R0R04	TI552 R0R04	TI752 R0R04
1000 mm	TI172 R0R04	TI372 R0R04	TI572 R0R04	TI772 R0R04

with double castors Ø 100 mm and drawer

- ▶ top shelf: flat, on tubular frame
- bottom shelf: flat, welded on



size ▶	650 x 450 mm	750 x 650 mm	950 x 650 mm	1250 x 650 mm
height ▼	Article number			
650 mm	TI112 F0154	TI312 F0154	TI512 F0154	TI712 F0154
800 mm	TI132 F0154	TI332 F0154	TI532 F0154	TI732 F0154
900 mm	TI152 F0154	TI352 F0154	TI552 F0154	TI752 F0154
1000 mm	TI172 F0154	TI372 F0154	TI572 F0154	TI772 F0154

with swivel castors Ø 125 mm and a push handle

- ▶ top shelf: with raised edge, on tubular frame
- bottom shelf: with raised edge, welded on





size ▶	650 x 450 mm	750 x 650 mm	950 x 650 mm	1250 x 650 mm
height ▼	Article number			
650 mm	TI112 R0327	TI312 R0327	TI512 R0327	TI712 R0327
800 mm	TI132 R0327	TI332 R0327	TI532 R0327	TI732 R0327
900 mm	TI152 R0327	TI352 R0327	TI552 R0327	TI752 R0327
1000 mm	TI172 R0327	TI372 R0327	TI572 R0327	TI772 R0327

INSTRUMENT TROLLEYS

Configure your own instrument table.



Select the required width and depth

T I

 width/depth	650 x 450 mm	750 x 650 mm	950 x 650 mm	1250 x 650 mm
	1	3	5	7



Select the required height

T I

 height	650 mm	800 mm	900 mm	1000 mm
	1	3	5	7

Select the number of shelves/tops



























T I

	2		3
--	---	--	---

top shelf middle shelf (0 with 2 shelves) bottom shelf

Select the type of shelf/top

T I

top				
		 1		 R
		 F		
center				
 0		 1		 3
		 F		 R
bottom				
		 1		 3
		 F		 R

Legend:



flat shelf (1 or F)

shelf with raised edge (3 or R)

shelf with gallery (G)



welded in place

removable

INSTRUMENT TROLLEYS

Configure your own instrument table.





Decide whether you need a drawer and/or a push handle

T I _____

0	2	5	7
without	with a push handle	with a drawer	with a push handle and a drawer
			

Choose suitable castors

T I _____

0	1	4	7
Ø 65 mm	Ø 75 mm	Ø 100 mm	Ø 125 mm
			

Enter the corresponding numbers and letters, one after another, into the respective fields and you will see the required article number.

T I _____

That's all you need to do! Simply enter the article number of your customized instrument trolley in your order.

Modular construction – 25,000 possible combinations – maximum freedom of choice.

Feel free to ask one of our persons of contact if you need help in configuring instrument or ISO trolleys to meet your special needs.

ISO TROLLEYS

QUALITY AND TECHNICAL DETAILS

OUR ARTICLE NUMBER: TM _ _ _ _ _

Tube

Module



42.

QUALITY

Made of chrom nickel steel tubes and sheets, material no. 1.4301 (AISI 304) with a 320 grit surface finish.

WASHABLE


ISO trolleys with these markings are suitable for machine washers.



CASTORS

All castors are of premium design: electrically conductive, non-marking and with thread guard.

lockable polymer double castors, Ø 100 mm



lockable, washable polymer swivel castors, Ø 125 mm



ISO TROLLEYS



Description

ISO trolleys 600 x 400 mm

- ▶ max. 6 trays 100 mm high
- ▶ with 6 layers
- ▶ top: with frame for holding a tray or chrome nickel steel plate with raised edge
- ▶ center: wire bracket for holding 4 trays
- ▶ bottom: shelf with raised edge, welded on
- ▶ with push handle

4 lockable polymer double castors
Ø 100 mm, wall protectors

Components

- ▶ Plate with raised edge for ISO-Table, chrome nickel steel

Article number

TM01064

W/D/H

870 x 535 x
1000 mm





















TM0180K

665 x 420 x
20 mm

Accessories TI und TM

	Description	Article number.:	 W/D/H
	<p>extension plate</p> <ul style="list-style-type: none"> ▶ additional option ▶ chrome nickel steel ▶ flat top ▶ can be folded down ▶ increases worktop area ▶ hygienic design 	TI901	350 x 650 x 30 mm
	<p>Standard rail structure for hanging TI and TM accessories on frame</p> <ul style="list-style-type: none"> ▶ additional option ▶ chrome nickel steel ▶ two standard rails ▶ additional adjustment and mounting options 	TI935	650 x 45 x 350 mm (the trolley is 45 mm wider and 380 mm higher)
 	<p>MAGNES side panel for mounting TI and TM accessories magnetically, can be removed without tools</p> <ul style="list-style-type: none"> ▶ MAGNES side panel ▶ MAGNES side panel with push-handle ▶ MAGNES side panel with standard rail 	<p>TI910</p> <p>TI912</p> <p>TI914</p>	<p>650 x 20 x 550 mm</p> <p>650 x 100 x 680 mm (the trolley is 80 mm wider and und 65 mm higher)</p> <p>650 x 70 x 550 mm (the trolley is 50 mm wider)</p>

Dispensers for syringes and cannulas

	Description	 W/D/H			
	3 tilting containers	600 x 215 x 245 mm	W1060 70010	W1460 70010A	
	4 tilting containers	600 x 180 x 210 mm	W1060 70100	W1460 70100A	
	5 tilting containers	600 x 145 x 165 mm	W1060 71000	W1460 71000A	
	5 + 3 tilting containers	600 x 200 x 420 mm		W1460 71010A	
	2 x 4 tilting containers	600 x 170 x 420 mm		W1460 70200A	
	5 + 4 tilting containers	600 x 170 x 390 mm		W1460 71100A	
	2 x 5 small tilting containers	600 x 170 x 340 mm		W1460 72000A	
	First in - First out	 W/D/H			
	3 FI/FO containers	600 x 220 x 245 mm	W1060 70011	W1460 70011A	
	5 FI/FO containers	600 x 165 x 165 mm	W1060 71001	W1460 71001A	
	5 + 3 FI/FO containers	600 x 210 x 420 mm		W1460 71011A	



Magnetic mounting on
MAGNES side panels









Products to be hung on the
standard rail attachment









Products to be hung on
a panel with standard rail

Accessories TI und TM

Glove dispensers

	Description	 W/D/H			
 RAL	Glove dispensers Hori 1 for 1 box of non-sterile gloves	max. box size 285 x 90 x 135 mm	W1030 20010	W1430 20010U	W1430 20010U
 RAL	Glove dispensers Vert 1 for 1 box of non-sterile gloves	max. box size 90 x 285 x 135 mm	W1030 20010	W1430 20015A	W1430 20015W






Disposal boxes

	Description	 W/D/H			
 CNS	Holder and set of 10 disposal boxes	à 1,5 l, Ø 140 x 145 mm, PP	W1050 40001	W1450 40001U	W1450 40001U
 CNS	Holder for box or container used by customer (when ordering, please state make and dimensi- ons of the boxes used)	up to 12 liters	W1050 40016	W1450 40016U	W1450 40016U








Disinfection

	Description	 W/D/H			
	Dispenser for one 500 ml Euro- format bottle of detergent/dis- infectant, bottle not included, manual operation		W1032 90050	W1432 90050A	W1432 90050W
	Dispenser for one 500 ml Euro- format bottle of detergent/dis- infectant, bottle not included, manual operation, with drip- catcher tray (CNS)		W1032 90052	W1432 90052A	W1432 90052W
 CNS	FlowPack holder for max. package size approx. 250 x 120 mm	120 mm high	W1038 40020	W1438 40020A	W1438 40020W





Waste bins

	Description	 W/D/H			
 RAL	Waste bin system with 19-litre bin (PP) and 1 roll of 35-litre waste bags	400 x 190 x 415 mm	W1055 20191		W1455 20191W

Disposal tray holder

	Description	 W/D/H			
 CNS	Holder for mounting PP container, mounting on long side of container		W1025 40014	W1425 40014U	W1425 40014U
 CNS	Holder for mounting a PP container, mounting on short side of container		W1025 40015	W1425 40015U	W1425 40015U
 CNS	Holder for mounting a container (GN 1/3)		W1025 40210	W1425 40210U	W1425 40210U

Disposal trays

	Description	 W/D/H	Article number.:
	Polypropylene container with slotted lid, sterilisable up to 134 °C	300 x 225 x 130 mm	W9150 70001
	Polypropylene container with lid, sterilisable up to 134 °C	300 x 225 x 130 mm	W9150 70002
 CNS	GN container / drip-catcher bowl GN 1/3-100, volume approx. 4 litres	325 x 176 x 100 mm	W9120 40210



Magnetic mounting on
MAGNES side panels



Products to be hung on the
standard rail attachment



Products to be hung on
a panel with standard rail

SERVICE TROLLEYS

QUALITY AND TECHNICAL DETAILS

OUR ARTICLE NUMBER: TG _ _ _ _ _

Tube

Gastronorm

48.



QUALITY

Made of chrom nickel steel tubes and sheets, material no. 1.4301 (AISI 304) with a 320 grit surface finish.

WASHABLE

Serving trolleys with these markings are suitable for machine washers.



CASTORS

All castors are of premium design: electrically conductive, non-marking and with thread guard.

lockable, washable polymer swivel castors, Ø 125 mm




BOARD/SHELVES


Deep-drawn raised edge, welded in place.

PUSH HANDLES

Two round tubular push handles Ø 25 mm.



H m. Rollen 125	1050	1250
h_1	200	100
h_2	565	
h_3	270	
h_4		270



B, T	400	550	750	1050
b, t	360	510	710	1010
GB, GT	490	640	840	1140

SERVICE TROLLEYS



Description

for 2 GN 1/1
NM 400 x 550

- ▶ 2 shelves with raised edge, welded in place

4 lockable, washable polymer castors Ø 125 mm, wall protectors

Article number.:

TG132 30307

W/D/H

490 x 640 x 1050 mm



Description

for 3 GN 1/1
NM 400 x 550

- ▶ 3 shelves with raised edge, welded in place

4 lockable, washable polymer castors Ø 125 mm, wall protectors

Article number.:

TG133 33307

W/D/H

490 x 640 x 1050 mm



Description

for 4 GN 1/1
NM 400 x 550

- ▶ 4 shelves with raised edge, welded in place

4 lockable, washable polymer castors Ø 125 mm, wall protectors

Article number.:

TG154 33307

W/D/H

490 x 640 x 1250 mm



Description

für 4 GN 1/1
NM 750 x 550

- ▶ 2 shelves with raised edge, welded in place

4 lockable, washable polymer castors Ø 125 mm, wall protectors

Article number.:

TG332 30307

W/D/H

840 x 640 x 1050 m



Description

für 6 GN 1/1
NM 750 x 550

- ▶ 3 shelves with raised edge, welded in place

4 lockable, washable polymer castors Ø 125 mm, wall protectors

Article number.:

TG333 33307

W/D/H

840 x 640 x 1050 mm

SERVICE TROLLEYS



Description

for 8 GN 1/1
NM 750 x 550

- ▶ 4 shelves with raised edge, welded in place

4 lockable, washable polymer castors Ø 125 mm, wall protectors

Article number.:

TG354 33307

W/D/H

840 x 640 x
1250 mm



Description

for 6 GN 1/1
NM 1050 x 550

- ▶ 2 shelves with raised edge, welded in place

4 lockable, washable polymer castors Ø 125 mm, wall protectors

Article number.:

TG532 30307

W/D/H

1140 x 640 x
1050 mm



Description

for 9 GN 1/1
NM 1050 x 550

- ▶ 3 shelves with raised edge, welded in place

4 lockable, washable polymer castors Ø 125 mm, wall protectors

Article number.:

TG533 33307

W/D/H

1140 x 640 x
1050 mm



Description

for 12 GN 1/1
NM 1050 x 550

- ▶ 4 shelves with raised edge, welded in place

4 lockable, washable polymer castors Ø 125 mm, wall protectors

Article number.:

TG554 33307

W/D/H

1140 x 640 x
1250 mm

GN = Gastronorm
NM = Nominal size

BECKER

WE SHINE WITH STEEL



**Metallwarenfabrik
Walter H. Becker GmbH**
Anzenkirchener Straße 4
84371 Triftern, Germany

Phone number +49 (0) 8562 9605-0
info@becker-triftern.de
www.becker-triftern.de

