

EES Skin Staplers

“I get versatility and a great value.”



Product Overview

The **PROXIMATE[®] RH** skin stapler is a full feature skin stapler with a 360° rotating head assembly. It features improved visibility over existing fixed head skin stapling products along with a ratchet mechanism for easier staple placement.

The **PROXIMATE[®] PX** Skin Stapler provides many of the same features as the PROXIMATE[®] RH skin stapler, but in a fixed head format.


The **PROXIMATE[®] PLUS MD** is a high-value, low-cost skin stapler that permits multi-directional release in an ergonomic design.

Product	Features and Benefits
PROXIMATE[®] RH	Stronger rectangular staples resist opening and minimize staple rotation, providing better cosmetic results Head rotates 360° for improved visibility and access Ergonomic pistol-grip handle is intuitive and comfortable to use Positive ratchet mechanism provides easy staple placement Staples are coated with Dry Film for easy staple extraction
PROXIMATE[®] PX	Ergonomic pistol-grip handle is intuitive and comfortable to use Positive ratchet mechanism provides easy staple placement Staples are coated with Dry Film for easy staple extraction
PROXIMATE[®] PLUS MD	Improved kick-off spring design for multi-directional release Ergonomic design comfortable for smaller hands Alignment indicator improves visibility Staples coated with Dry Film for easy, less painful staple extraction




Ethicon
Endo-Surgery


PROXIMATE® RH

Rotating Head Skin Stapler	Product Code	Type	No of Staples		per pkg.	Staple Dimensions (mm)		
						Wire Diam.	Leg Length	Crown Length
	PRR35	Regular	35		6	0.53	3.9	5.7
	PRW35	Wide	35		6	0.58	3.9	6.9

PROXIMATE® PLUS MD

Multi-Directional Release Skin Stapler	Product Code	Type	No of Staples		per pkg.	Staple Dimensions (mm)		
						Wire Diam.	Leg Length	Crown Length
	PMR35	Regular	35		6	0.53	3.9	5.7
	PMW35	Wide	35		6	0.58	3.9	6.9

PROXIMATE® PX

Skin Stapler	Product Code	Type	No of Staples		per pkg.	Staple Dimensions (mm)		
						Wire Diam.	Leg Length	Crown Length
	PXR35	Regular	35		6	0.53	3.9	5.7
	PXW35	Wide	35		6	0.58	3.9	6.9

Miscellaneous Products	PSX	Squeeze-Handle Staple Extractor (for all models)						
------------------------	-----	--	--	--	--	--	--	--

